


RADIOFREQUENCY ABLATION OF THYROID NODULES

Procedure Guide



DISCLAIMER

The surgical technique herein is presented to demonstrate the method utilized by Prof. Dr. Dr. Hüdayi Korkusuz. The information on the products and procedures contained in this brochure does not represent and does not constitute medical advice or recommendations and should not be relied upon as such.

This information does not purport to constitute any diagnostic or therapeutic statement with regard to any individual medical case. Each patient must be examined and advised individually, and this brochure does not replace the need for such examination and/or advice in whole or in part.

This brochure should not be considered as a substitute for carefully reading all applicable labeling, including the instructions for use (IFU) supplied with the devices. Before using any product, please thoroughly review the relevant user manual(s) for instructions, including, but not limited to, contraindications, warnings, precautions, and adverse effects. Please note: It is the clinicians' responsibility to decide which instrument mode and settings they use in each clinical situation.

- Bipolar RFA is a safe and effective treatment option for symptomatic benign thyroid nodules.(1)
- Ultrasound-guided percutaneous BRFA seems to be an effective and safe method for the treatment of benign thyroid nodules. It may gain a wide use in clinical practice.(2)
- The use of bipolar RFA is an effective, safe and suitable thermoablative technique to treat benign thyroid nodules. Combined with the multiple overlapping shot technique it allows sufficient ablation.(3)

PREPARATION AND ANESTHESIA

Radiofrequency ablation of the thyroid nodule

In this application guide the treatment of benign thyroid nodules with the aid of bipolar radiofrequency ablation is demonstrated. Before you use this method, make sure that a proper diagnosis of the type of the nodule has been made and that malignancy has been excluded e.g. by fine needle aspiration.

1



The patient shall lie down comfortably on a treatment table with the neck extended slightly.

6

- Before starting the procedure, the location of the thyroid nodule is confirmed with the aid of sonographic imaging.
- The dimensions of the nodule are measured in transisthmic view in order to choose a suitable applicator for the treatment.

Example:

Nodule diameter: 32 mm

Chosen applicator: CELON ProSurge 100-T20 (20 mm electrode length)

Preparation and Anesthesia

4 After choosing an appropriate applicator, the generator can be started and the CELON ProSurge applicator is connected. More information about how to prepare the system can be found in our "Olympus CELON Power System Set-Up Guide".



5


disinfection to the patient's skin in order to avoid infections.

Under ultrasound guidance, local

It is important

to provide an

extensive amount of




of the thyroid or the thyroid nodule is generally not necessary. In some cases, it might be advisable to administer an analgesic to your patient. You can use a scalpel to make a 2 mm incision in the skin to ease the insertion of the CELON ProSurge applicator.

For the insertion of the applicator, sterile conditions are mandatory.

TREATMENT OF THE THYROID

9





By applying the transisthmic approach, vulnerable structures can be displayed and avoided.

Under ultrasound guidance, the applicator is now inserted from the contralateral side through the isthmus of the thyroid into the nodule.

Treatment of the Thyroid


A safety margin next to the so called danger triangle containing the recurrent laryngeal nerve between trachea and thyroid should be created.

You should also provide a safety margin to vulnerable structures like the carotid artery, the jugular vein and the vagus nerve.

Internal jugular vein

Vagus nerve (next to internal carotid artery)

Common carotid artery



Thyroid nodule

Danger Triangle incl. recurrent laryngeal nerve


Illustration of vulnerable structures alongside a thyroid nodule being visible under ultrasound

By inserting the applicator in the transisthmic approach, the entire length of the electrode can be visualized via ultrasound.

The insulator between the two electrodes at the tip is clearly visible.

TREATMENT OF THE THYROID

- After positioning the applicator, activate the RF power by pressing the foot switch. You can monitor the status of coagulation via sonographic imaging.
- On the ultrasound display, an upcoming hyperechogenic cloud becomes visible around the insulator. It can be identified as an indicator for vaporizing.



Hyperechogenic cloud

CELON ProSurge applicator

16

Ultrasound image of a CELON ProSurge applicator inside a thyroid nodule.

The ablation of the first part of the thyroid nodule is finished when the power output slightly decreases.

The generator power display switches from set power to applied power after pressing the foot switch.

Electrode	Set power	Target power
T20	25W	15W
.20		
T30	35W	25W
T40	45W	35W

It is recommended to animate the patient to a conversation during the ablation to recognize a possible damage of the nerves in an early stage.

Treatment of the Thyroid

In the following step, the position of the applicator is slightly changed and the ablation is started again.

By repositioning the applicator next to the prior ablation area, you can create overlapping ablation volumes. With only a few repositionings, most of the nodule volume is treated. This approach is called Multiple Overlapping Shot Technique (MOST).



- This step is repeated until the required ablation volume is achieved. An indicator for complete treatment is total coverage of the targeted nodule volume by the hyperechogenic cloud.
- 19 After the ablation procedure is completed, the applicator can be easily removed by pulling it out gently.
- The insertion site just needs to be cleaned with a swab.

THERAPEUTIC EFFECT

21



The coagulation achieves a local denaturation of the treated area.

Therapeutic Effect

- 22 After several weeks, your patient can expect a visible reduction in volume as a result of the body's resorption of the coagulated tissue.
- The therapeutic effect can also be measured with the help of different imaging modalities like echogenicity and elasticity measurements, contrast-enhanced ultrasound, Doppler blood flow and (MIBI) scintigraphy.



POSTOPERATIVE CARE

24



Immediately after the treatment, the patient will hardly notice any differences. The puncture site only needs to be covered with a band-aid.

Postoperative Care

- Following the procedure, the patient should help protect the wound.

 Advise your patient not to take showers over the next couple of days.

 Furthermore, overstretching of the neck as well as heavy physical work should be avoided.
- Schedule a follow-up appointment with your patient after three months to control the effect of the bipolar radiofrequency ablation of the thyroid nodule.

RADIOFREQUENCY ABLATION OF THYROID NODULES

Ordering Information

CELON Power System

WB992001 Set "CelonPOWER System"

Consisting of generator CelonLab POWER (WB991029), peristaltic pump CelonAquaflow III (WB950059) and system

trolley CelonMobile (WB950067)

CELON ProSurge Applicators

Internally cooled bipolar radiofrequency applicator for tumor ablation to be used individually or in multipolar combinations.

Included in delivery of each applicator: Applicator, tubing, puncture guide

1.8mm (15 Gauge), three-sided tip, sterile, for single use

Order Number	Description	Shaft lengths/electrode lengths
WB990129	CelonProSurge 100-T20	100 mm / 20 mm
WB990146	CelonProSurge 100-T30	100 mm / 30 mm
WB990147	CelonProSurge 100-T40	100 mm / 40 mm
WB990148	CelonProSurge 150-T20	150 mm / 20 mm
WB990149	CelonProSurge 150-T30	150 mm / 30 mm
WB990150	CelonProSurge 150-T40	150 mm / 40 mm
WB990151	CelonProSurge 200-T20	200 mm / 20 mm
WB990152	CelonProSurge 200-T30	200 mm / 30 mm
WB990153	CelonProSurge 200-T40	200 mm / 40 mm
WB990187	CelonProSurge 250-T30	250 mm / 30 mm
WB990188	CelonProSurge 250-T40	250 mm / 40 mm



CELON ProSurge micro Applicators

Bipolar RF applicator for precise thermal ablation of tissue

1.3mm (18 Gauge), three-sided tip, sterile, for single use

Order Number	Description	Shaft lengths/electrode lengths
WB990072	CelonProSurge micro 100-T09	100 mm / 9 mm
WB990104	CelonProSurge micro 100-T15	100 mm / 15 mm
WB990091	CelonProSurge micro 150-T09	150 mm / 9 mm
WB990105	CelonProSurge micro 150-T15	150 mm / 15 mm
WB990092	CelonProSurge micro 200-T09	200 mm / 9 mm
WB990106	CelonProSurge micro 200-T15	200 mm / 15 mm



References

- 1 Korkusuz Y, Erbelding C, Kohlhase K, et al. Bipolar Radiofrequency Ablation of Benign Symptomatic Thyroid Nodules: Initial experience with Bipolar Radiofrequency. Fortschr Röntgenstr 2016; 188: 671–675.
- 2 Li XL, Xu HX, Lu F, et al. Treatment efficacy and safety of ultrasound-guided percutaneous bipolar radiofrequency ablation for benign thyroid nodules; Br J Radiol. 2016 Mar; 89(1059):20150858.
- 3 Kohlhase KD, Korkusuz Y, Gröner D, et al. Bipolar radiofrequency ablation of benign thyroid nodules using a multiple overlapping shot technique in a 3-month follow-up. Int J Hyperthermia. 2016 Aug;32(5):511-6.

