


MECHANICALLY ENHANCED COLONOSCOPY  
**ENDOCUFF VISION™**  
Master the Art of Unfolding




# Master the Art of Unfolding

with ENDOCUFF VISION™

## The Colon Is Complex

Tortuous anatomy can turn a routine colonoscopy into a series of challenges.



## Increase Your ADR


Designed to increase adenoma detection rates (ADR), ENDOCUFF VISION™ maximizes visualization by flattening mucosal folds while helping improve control during colonoscopy.

## Detection Means Prevention

Increasing the Adenoma Detection Rate Improves Your Ability to Save Lives

“ ADR is considered the primary measure of the quality of mucosal inspection and the single most important quality measure in colonoscopy. ”

American Society for Gastrointestinal Endoscopy, Quality Indicators for Colonoscopy 2015<sup>9</sup>



- A 1% increase in ADR results in a 3% decrease in the risk of interval cancer and a 5% decrease in the risk of a fatal interval colorectal cancer<sup>5</sup>
- Increasing ADR not only lowers mortality rates but also reduces treatment costs<sup>6</sup>

### Unseen Adenomas and Polyps Are Causes for Concern




**24%** of adenomas are missed with standard colonoscopy<sup>7</sup>



**40%** As many as 40% of polyps go undetected<sup>8</sup>

## Vision Leads to Action

Providing a Better View of the Colon by Flattening Mucosal Folds



#### Control:

- Anchors endoscope tip in the lumen during examination and polypectomy.
- Reduces slippage and minimizes difficulties associated with looping.

#### Application:

- Hard plastic body and firm, dry grip on the scope prevent dislodgement and protect distal tip of the colonoscope.
- Hinged arms fit seamlessly into device to create a low-profile design.

#### Design:

- The proprietary hinged design allows the arms to fall flat against the scope for smooth forward movement during intubation.

#### Visualization:

- Upon withdrawal, the uniquely hinged arms expand to gently flatten large mucosal folds, bringing difficult to see areas into view.
- Soft, flexible arms provide the right amount of force to be effective without causing mucosal trauma.

# Obtain Visible Results


## Clinical Studies Show Greater Visibility

### Screening Performance Shows Greater Visibility Compared with Standard Colonoscopy in Evaluation Studies

- The use of ENDOCUFF VISION™ resulted in a notable increase in ADR of 10.8%.<sup>10</sup>
- ENDOCUFF VISION™ appeared to improve operator performance without more patient discomfort through shorter procedure times and no requirement for higher sedation.<sup>7</sup>



Adenoma Detection Comparison (p<0.001)



### Lesion Detection Compared with Standard Colonoscopy in the BCSP Population

	Standard Colonoscopy (n=403)	ENDOCUFF VISION™ Assisted Colonoscopy (n=394)	p-Value
Mean Adenomas per Procedure	1.20	1.59	p=0.004
Polyps	63.3%	73.9%	p<0.001
Sessile Serrated Adenomas	1.2%	2.0%	p=0.19
Left Colon Adenomas	32.8%	40.9%	p=0.009
Right Colon Adenomas	38.0%	43.2%	p=0.07
Large Adenomas (>10 mm)	12.4%	13.7%	p=0.29
Small Adenomas (6-9 mm)	10.7%	19.0%	p<0.001
Diminutive Adenomas (≤5 mm)	44.7%	52.0%	p=0.02

Adapted from Ngu et al. 2018.<sup>11</sup>



During intubation, the flexible arms of ENDOCUFF VISION™ slip into the body of the device so that forward movement is not hindered.



The arms open folded bowel for inspection, everting large mucosal folds and providing a clear view of the mucosa that was previously difficult to visualize.



During withdrawal, the flexible arms reduce the risk of sudden slippage, stabilize the scope tip during therapy, and anchor the scope tip during loop reduction.

### Order Information

Article Name	Article Number	Description	Compatible Scope(s)	Quantity Units/Box
ARV110	K10037616	ENDOCUFF VISION™ Medium Blue I.D. 11.0	PCF-H190DL/I, CF-H185L/I, CF-Q145L/I, CF-H290L/I, CF-Q260AL/I, CF-Q260DL/I, CF-240AL/I, CF-240DL/I, CF-240L/I	8
ARV120	K10037617	ENDOCUFF VISION™ Large Green I.D. 11.2	CF-HQ190L/I, CF-H190L/I, CF-H180AL/I, CF-Q180AL/I, CF-Q165L/I, CF-Q160AL/I, CF-Q160DL/I, CF-Q160L/I, CF-Q160S, CF-Q140L/I, CF-140L/I, CF-140S, CF-1T140L/I, CF-HQ290L/I, CF-H260AL/I, CF-H260DL/I, CF-H260AZL/I, CF-Q240AL/I, CF-Q240L/I, CF-Q240ZL/I, CF-H170L/I	8
ARV130	K10037618	ENDOCUFF VISION™ Small Purple I.D. 10.4	PCF-H190L/I, PCF-H180AL/I, PCF-Q180AL/I, PCF-160AL/I, PCF-H290L/I, PCF-Q260AL/I, PCF-Q260AZI	8
ARV140	K10037619	ENDOCUFF VISION™ XL Orange I.D. 12.1	CF-H180DL/I	8

[www.olympus.eu/et-catalog](http://www.olympus.eu/et-catalog)

### References:

- Witte TN, Enns R. The difficult colonoscopy. *Can J Gastroenterol.* 2007;21(8):487-490.
- Rex DK, Ahnen DJ, Baron JA, et al. Serrated Lesions of the Colorectum: Review and Recommendations From an Expert Panel. *Am J Gastroenterol.* 2012;107(9):1315-1329.
- Pickhardt PJ, Nugent PA, Mysliwiec PA, et al. Location of Adenomas Missed by Optical Colonoscopy. *Ann Intern Med.* 2004;141(5):352-359.
- Jayasekaran V, Holt B, Bourke M. Normal Adult Colonic Anatomy in Colonoscopy. *Video Journal and Encyclopedia of GI Endoscopy.* 2013;1(2):390-392.
- Corley DA, Jensen CD, Marks AR, et al. Adenoma Detection Rate and Risk of Colorectal Cancer and Death. *N Engl J Med.* 2014;370(14):1298-1306.
- Meester RG, Doubeni CA, Lansdorp-Vogelaar I, et al. Variation in Adenoma Detection Rate and the Lifetime Benefits and Cost of Colorectal Cancer Screening: A Microsimulation Model. *JAMA.* 2015;313(23):2349-2358.
- Tsiamoulos ZP, Ravi M, et al. Impact of a new distal attachment on colonoscopy performance in an academic screening center. *GIE Journal.* 2017.
- Steele SR, Johnson EK, Champagne B, et al. Endoscopy and polyps — diagnostic and therapeutic advances in management. *World J Gastroenterol.* 2013;19(27):4277-4288.
- Rex DK, Schoenfeld PS, Cohen J, et al. Quality indicators for colonoscopy. *Gastrointest Endosc.* 2015;81(1):31-53.
- Floer M, Biecker E, Fitzlaff R, et al. Higher Adenoma Detection Rates with Endocuff-Assisted Colonoscopy — A Randomized Controlled Multicenter Trial. *PLoS ONE.* 2014;9(12):e114267.
- Ngu WS, Bevan R, Tsiamoulos ZP, et al. Improved Adenoma Detection with Endocuff Vision: The ADENOMA randomized controlled trial. *GUT* 2018;66:1-9.

# ENDOCUFF VISION™

Olympus is a registered trademark of Olympus Corporation, Olympus America Inc. and/or their affiliates.

As medical knowledge is constantly growing, technical modifications or changes of the product design, product specifications, accessories and service offerings may be required.

**OLYMPUS**

**OLYMPUS EUROPA SE & CO. KG**

Postbox 10 49 08, 20034 Hamburg, Germany  
Wendenstrasse 20, 20097 Hamburg, Germany  
Phone: +49 40 23773-0, Fax: +49 40 233765  
[www.olympus-europa.com](http://www.olympus-europa.com)