

Spiration Valve System

A Proven Concept in Severe Emphysema Treatment



Severe Emphysema

A Progressive Disease with a Groundbreaking Treatment Horizon

Emphysema is a type of chronic obstructive pulmonary disease (COPD) that is progressive in nature and characterized by loss of elasticity and enlargement of the alveolar space in the lung. As a result, the diseased portion of the lung becomes hyperinflated, causing significant breathing challenges.



Symptoms of Severe Emphysema

- Breathlessness
- Fatigue
- Limitations to daily activities
- Reduced quality of life
- Reduced life expectancy

Current Treatment Options

- Smoking cessation
- Medical management
- Pulmonary rehabilitation
- Oxygen therapy
- Surgical intervention



Based on safety and efficacy data from multiple international clinical studies- bronchoscopic lung volume reduction (BLVR), using endobronchial valves (EBV), is now recommended by numerous prominent guidelines as a treatment option for advanced emphysema.

These include:

- Disease (GOLD)¹
- National Institute for Health and Care Excellence Interventional Procedures Guidance (NICE)²

Endobronchial Valve Treatment

Now a globally recommended therapy

2019 Global Initiative for Chronic Obstructive Lung

The Spiration Valve System

The Right Patient. The Right Valve. The Right Outcomes.

The Spiration Valve System (SVS) for bronchoscopic lung volume reduction is proven to improve lung function, reduce shortness of breath, and restore quality of life.³

The SVS has demonstrated a strong risk benefit profile, with a low rate of serious pneumothorax and minimal risk of valve migration and expectoration.*

It also offers noninvasive patient selection, a short procedure time, and the assurance of Olympus' expertise in medical and respiratory technology.



 * Please refer to full prescriptive information in the back of this document



The Spiration Valve System

Innovative Technology

The Spiration Valve System (SVS) is an innovative endobronchial technology that offers patients with severe emphysema a customized, minimally-invasive treatment option for lung volume reduction with a favorable risk-benefit profile.



The Spiration Valve is delivered to the target lobe during a bronchoscopic procedure. The valve anchors are designed to maintain position and minimize expectoration.



On inhalation the valve redirects air to healthier **portions** of the lung enabling healthier tissue to expand.



On exhalation the valve allows trapped air and secretions to escape from the hyperinflated lobe.

By allowing air to leave but not re-enter diseased areas of the lung, it is possible to reduce hyperinflation in the targeted lobe. Treatment typically requires placement of multiple valves to achieve complete lobar occlusion in targeted lobe.

Improve lung function

Bronchoscopic Lung Volume Reduction With the SVS Has Been Shown To:

- Allow healthier lung to re-expand
- Reduce dyspnea Increase the ability to perform daily activities



In clinical trials, patients treated with the SVS experienced improved breathing, lung function, and quality of life.³







Procedure Overview

A short bronchoscopic procedure under general anesthesia or deep sedation.

Airway Sizing and Valve **Selection for a Custom Fit**

A calibrated balloon is used to customize the appropriate valve size for placement.



Easy, Reliable Valve Loading and Deployment

The cartridge and loader work seamlessly to quickly compress the valve into the catheter in preparation for valve placement.



Versatile Valve Placement

Flexible valve design and delivery catheters enable placement in targeted airways.



Easy Valve Removal

Valve designed to facilitate retrieval with 360° access to removal rod.

The Spiration Valve

Engineered for Dynamic Lung Anatomy



escape naturally along the

bronchial wall.



A center rod of the valve facilitates easy removal.

Design Overview

The Spiration Valve is designed to allow flexible placement even in tortuous anatomy, such as the airways in the upper lobe segments of the lungs.



Multiple valve sizes accommodate variable lung anatomy with precision fit.



Evidence Confirms Effectiveness³

EMPROVE Clinical Trial

Study Summary⁴

Randomized, controlled clinical trials (RCTs) demonstrate the SVS is a safe and effective treatment option for patients with severe emphysema and little-to-no collateral ventilation.



Study Overview

- with severe emphysema
- 2:1 randomization into SVS treatment arm (n=113), and standard of care control arm (n=59) Alpha-1 antitrypsin deficiency non-randomized SVS treatment arm (n=20) Primary and secondary effectiveness endpoints measured six months following randomization

- Longer term durability of effectiveness measured at 12 months following randomization

Patient Selection

EMPROVE exclusively used high-resolution computed tomography (HRCT), a non-invasive approach to identifying patients with low to no collateral ventilation.

Summary of Outcomes^{3,4}

The EMPROVE clinical trial demonstrated that patients treated with the SVS benefited from significant clinical and statistical improvements in lung function and quality of life, compared to standard of care medical management.

PRIMARY ENDPOINT

Change in FEV,



* A negative change in SGRQ represents an improvement in disease specific health status. A 4 point reduction is considered clinically meaningful.

EMPROVE evaluated the safety and effectiveness of the Spiration Valve System in 172 patients

Evidence Confirms Effectiveness³

EMPROVE Clinical Trial

SECONDARY ENDPOINTS

Targeted Lobe Volume Reduction³



Pre-Treatment: 1.84 <u>+</u> 0.6 L



53%

PP: 1.0000

Reduction

Difference: -0.974 L

95% BCI: -1.119, -0.829

Post-Treatment: 0.87 <u>+</u> 0.9 L

Atelectasis		
Complete (≥99% TLVR)	40/102 (39%)	
Substantial (≥50% TLVR)	51/102 (50%)	
MCID Responder (≥350ml TLVR)	76/102 (75%)	

Pulmonary Adverse Events ³	SVS Group % (N = 113) Short-Term	Control Group % (N = 59) (0-6) Months	SVS Group % (N = 113) Long-Term	Control Group % (N = 47) (6-12) Months
Acute exacerbation of COPD	16.8	10.2	13.6	8.5
Death from procedure or device	0.0	_	1.0	_
Pneumonia - in the valve-treated lobe	1.8	_	1.0	_
Pneumonia - in the non-valve-treated lobe	7.1	1.7	7.8	2.1
Serious Pneumothorax	14.2	0.0	0.0	0.0
Respiratory failure	2.7	0.0	1.0	0.0

St. George's Respiratory Questionnaire³



Dyspnea Score (mMRC)³



Hyperinflation (RV/TLC)³



The SVS treatment group showed statistically significant improvement at both 6 and 12 months compared to the control group. A 4 point reduction is considered clinically meaningful.

Mean ±95% Bayesian Credible Interval PP = Posterior Probability

Patient Selection

The Key to Successful Outcomes

A decade of clinical studies shows appropriate patient selection to be one of the most important predictive factors of an effective response to bronchoscopic lung volume reduction.⁶

A thorough patient evaluation, examination for any comorbidities, and analysis of the patient's HRCT information and quantitative computed tomography (QCT) results are critical to successful outcomes. The below criteria may be used as a guide for appropriate patient selection based on the EMPROVE Trial.*

EMPROVE Inclusion Criteria³

Assessment	Inclusion Criteria		
Medical History and Physical Exam	 ≥ 40 years of age Diagnosed with severe emphysema Considered to have "stable" COPD as defined by Guidelines for Management of Stable COPD⁵ ≥ 6 weeks without exacerbation Able to tolerate a bronchoscopic procedure 		
Radiographic	Severe emphysema defined as target lobe with $\ge 40\%$ emphysema involvement High heterogeneity defined as ≥ 10 point disease severity difference with the ipsilateral lobe Fissure integrity defined as $\ge 90\%$ completeness of the fissure(s) separating the target lobe		
Pulmonary and Exercise Evaluation	FEV ₁ Residual Volume (RV) Total Lung Capacity (TLC) 6MWD	< 45% predicted ≥150% predicted ≥100% predicted ≥140 meters	

Exclusion Criteria

- Patient is an active smoker.
- Patient has a severe gas exchange abnormality in either PCO, or PO, as defined by PCO, >55 mm Hg, or PO, < 45 mm Hg on room air.
- Patient has a BMI < 15kg/m².
- Patient had a hospitalization for COPD exacerbation or respiratory infections in the past 3 months prior to baseline testing.
- Patient has bronchitis with sputum production > 4 tablespoons per 60 ml per day.
- Patient has an active asthma component to their disease or requires more than 15mg of prednisone daily.

- Patient has giant bulla considered to be > 1/3 volume in either lung.
- Patient has severe pulmonary hypertension based upon clinical evaluation.
- Patient has had prior lung volume reduction surgery or major lung procedures (lobectomy or greater).
- Patient has a diffuse emphysema pattern.
- Patient is classified as ASA Class greater than P4 including presence of co-morbidity that could significantly increase the risk of a bronchoscopy procedure.³





* These recommendations are not meant to replace patient-specific clinical judgement.

14

SeleCT

SeleCT is a completely noninvasive patient screening solution that provides key measures of emphysema severity, fissure integrity, and heterogeneity.

These measures are provided in an easy-to-read report to assist with selecting qualified patients and potential target lobes for improved outcomes using bronchoscopic lung volume reduction (BLVR).

Key quantitative measures to identify responders for the Spiration Valve System:

M EMPHYSEMA SEVERITY

Allows physician to quickly identify the most diseased lobe

HETEROGENEITY

Differentiates target and ipsilateral lobe emphysema to facilitate redirection of ventilation to healthier tissue^{6,7}

FISSURE INTEGRITY

EMPROVE trial results confirmed radiographic assessment of fissure completeness to be a reliable surrogate for collateral ventilation³

The Most Comprehensive Solution

For Minimally Invasive Bronchoscopic Lung Volume Reduction

The Olympus Solution offers not only the state-of-the-art valve technology you need for effective BLVR, but also provides an entire portfolio of respiratory devices and bronchoscopes that ensure improved efficiency and quality patient care in the bronchoscopy suite.



EVIS X1 **Imaging Platform**

EVIS X1 introduces a range of new, easy-to-use technologies that combines high imaging performance and ergonomic working in daily routine.





The Spiration Valve is an umbrella-shaped, one-way valve that redirects air away from diseased area of the lung to healthier tissue, all while allowing trapped air and secretions to escape, so that patients can breathe easier.

Spiration Valve System Accessories

SVS offers a complete selection of supporting devices to ensure that valves can effectively be placed into the target airway, precisely deployed, and easily removed whenever deemed necessary.

SeleCT Quantitative CT Analysis

The SeleCT QCT service offers rapid results, including a qualified over-read by a certified thoracic radiologist.



BF-1TH1100 Bronchoscope

Fully rotatable, therapeutic bronchoscope with superb HDTV image quality and a 3.0mm working channel. This enables increased suction capabilities and contributes to better view in the bronchial airways.



Olympus Respiratory Solutions

Olympus offers a versatile product line-up from EndoTherapy devices and bronchoscopes to reprocessing and service solutions. Olympus is a trusted partner in the bronchoscopy suite when it comes to endoscopic interventions.











Spiration Valve System

References

- 2019 Global Strategy for the Diagnosis, Management, and Prevention of COPD. Global Initiative for Chronic Obstructive Lung Disease (GOLD) 2018. http://goldcopd.org
- National Institute for Health and Care Excellence (2017) Endobroncial valve insertion to reduce lung volume in emphysema (NICE guideline IPG 600) Available at: nice.org.uk/guidance/ipg600 [Accessed 16 November 2018].
- Criner GJ, Delage A, Voelker K. Late Breaking Abstract Endobronchial Valves for Severe Emphysema – 12-month Results of the EMPROVE Trial. Eur Respir J. 2018;52 (suppl 62). doi:10.1183/13993003.congress-2018.OA4928
- 4. Spiration Valve System. Summary of Safety and Effectiveness; 2018. Data on File.
- ACP/ACCP/ATS/ERS Guidelines for Management of Stable COPD. http://www.thoracic.org/statements/ resources/respiratory-disease-adults/ 179full.pdf. Accessed 2018.
- Eberhardt R, Gompelmann D, Schuhmann M, Heussel CP, Herth FJF. Complete unilateral vs partial bilateral endoscopic lung volume reduction in patients with bilateral lung emphysema. Chest. 2012;142(4):900-908. doi:10.1378 chest.11-2886
- 7. Thurnheer R, Engel H, Weder W, et al. Role of lung perfusion scintigraphy in relation

to chest computed tomography and pulmonary function in the evaluation of candidates for lung volume reduction surgery. Am J Respir Crit Care Med 1999;159(1):301-10.

 Wood DE, Nader DA, Springmeyer SC, et al. The IBV Valve trial: a multicenter, randomized, double-blind trial of endobronchial therapy for severe emphysema. J Bronchology Interv Pulmonol. 2014;21(4):288-297. doi:10.1097/ LBR.00000000000110



As medical knowledge is constantly growing, technical modifications or changes of the product design, product specifications, accessories and service offerings may be required.



OLYMPUS EUROPA SE & CO. KG Postbox 10 49 08, 20034 Hamburg, Germany Wendenstrasse 20, 20097 Hamburg, Germany Phone: +49 40 233773-0, Fax: +49 40 233765 www.olympus-europa.com