

## ENDOCUFF VISION™ Compatibility Chart\*

### OLYMPUS

| Series               | Scope Model   | ENDOCUFF VISION™ Model |
|----------------------|---------------|------------------------|
| EVIS X1              | CF-EZ1500DL/I | ARV120 (Green)         |
|                      | CF-XZ1200L/I  | ARV120 (Green)         |
|                      | CF-HQ1100DL/I | ARV120 (Green)         |
| EVIS EXERA III       | CF-HQ190L/I   | ARV120 (Green)         |
|                      | CF-H190L/I    | ARV120 (Green)         |
|                      | PCF-H190DL/I  | ARV110 (Blue)          |
|                      | PCF-H190L/I   | ARV130 (Purple)        |
|                      | PCF-HQ190L/I  | ARV110 (Blue)          |
|                      | PCF-H190TL/I  | Not compatible         |
|                      | CF-H185L/I    | ARV110 (Blue)          |
| EVIS EXERA II        | CF-H180AL/I   | ARV120 (Green)         |
|                      | CF-H180DL/I   | ARV140 (Orange)        |
|                      | CF-Q180AL/I   | ARV120 (Green)         |
|                      | PCF-Q180AL/I  | ARV130 (Purple)        |
|                      | PCF-H180AL/I  | ARV110 (Blue)          |
|                      | CF-Q165L/I    | ARV120 (Green)         |
| EVIS EXERA           | CF-Q160AL/I   | ARV120 (Green)         |
|                      | CF-Q160DL/I   | ARV120 (Green)         |
|                      | CF-Q160L/I    | ARV120 (Green)         |
|                      | CF-Q160S      | ARV120 (Green)         |
|                      | PCF-160AL/I   | ARV130 (Purple)        |
|                      | CF-Q145L/I    | ARV110 (Blue)          |
| EVIS 140             | CF-Q140L/I    | ARV120 (Green)         |
|                      | CF-140L/I     | ARV120 (Green)         |
|                      | CF-140S       | ARV120 (Green)         |
|                      | CF-1T140L/I   | ARV120 (Green)         |
|                      | PCF-140L/I    | Not compatible         |
| EVIS LUCERA ELITE    | CF-HQ290L/I   | ARV120 (Green)         |
|                      | CF-H290L/I    | ARV110 (Blue)          |
|                      | PCF-H290L/I   | ARV130 (Purple)        |
|                      | CF-H290ECI    | ARV120 (Green)         |
| EVIS LUCERA SPECTRUM | CF-H260AL/I   | ARV120 (Green)         |
|                      | CF-H260DL/I   | ARV120 (Green)         |
|                      | CF-H260AZL/I  | ARV120 (Green)         |
|                      | CF-Q260AL/I   | ARV110 (Blue)          |
|                      | CF-Q260DL/I   | ARV110 (Blue)          |
|                      | PCF-Q260AL/I  | ARV130 (Purple)        |
|                      | PCF-Q260AZI   | ARV130 (Purple)        |
|                      | PCF-PQ260L/I  | Not compatible         |
| EVIS LUCERA          | CF-Q240AL/I   | ARV120 (Green)         |
|                      | CF-Q240L/I    | ARV120 (Green)         |
|                      | CF-Q240ZL/I   | ARV120 (Green)         |
|                      | CF-240AL/I    | ARV110 (Blue)          |
|                      | CF-240DL/I    | ARV110 (Blue)          |
|                      | CF-240L/I     | ARV110 (Blue)          |
| OPTERA               | CF-H170L/I    | ARV120 (Green)         |

\* The table provides an indication of ENDOCUFF VISION™ compatibility with colonoscopes. For confirmation of specific models, contact an Olympus representative.

## ENDOCUFF VISION™ Compatibility Chart\*

### FUJI

| Scope Model                    | ENDOCUFF VISION™ Model |
|--------------------------------|------------------------|
| EC-760R-V/L Adult              | ARV110 (Blue)          |
| EC-760ZP-VM/VL Adult           | ARV110 (Blue)          |
| EC-600WM/WI/WL New Series      | ARV110 (Blue)          |
| EC-590WM4/WI4/WL4/WI/WLY Adult | ARV110 (Blue)          |
| EC-590ZW3 M/L Adult Zoom       | ARV120 (Green)         |
| EC-530FI/FL/FM Adult           | ARV120 (Green)         |
| EC-530WM3/WI3/WL3/WI           | ARV120 (Green)         |
| EC-450 Adult Series            | ARV120 (Green)         |
| EC-200MR Adult                 | ARV120 (Green)         |
| EC-250 Adult Series            | ARV120 (Green)         |
| EC-201WI Adult                 | ARV120 (Green)         |
| EC-530 MT/IT/LT Adult          | ARV120 (Green)         |
| EC-530 DM/DL Adult             | ARV120 (Green)         |
| EC-600 ZW/M/L Adult            | ARV120 (Green)         |
| EC-530MP/LP Pediatric          | ARV130 (Purple)        |
| EC-530LS Intermediate          | ARV130 (Purple)        |
| EC-450 Pediatric Series        | ARV130 (Purple)        |
| EC-450LP/LS5 Intermediate      | ARV130 (Purple)        |
| EC-250 Intermediate Series     | ARV130 (Purple)        |

### PENTAX

|                                       |                 |
|---------------------------------------|-----------------|
| EC38-i10M/M2/F/F2/L/NF Adult          | ARV110 (Blue)   |
| EC-3890TMK/MK2/FK/FK2/LK Adult        | ARV110 (Blue)   |
| EC-3890MZI/FZI/LZI/SI2/TFK Adult      | ARV110 (Blue)   |
| EC-3890FI/LI/Fi2(HD)/Mi2 Adult        | ARV110 (Blue)   |
| EC-3890TLK Dual Channel Adult         | ARV110 (Blue)   |
| EC-3970LK Adult                       | ARV120 (Green)  |
| EC-3880FK2 Adult                      | ARV120 (Green)  |
| EC-3870LK/LZK/TLK/FK2/MK2 Adult       | ARV120 (Green)  |
| EC-3840FK2 Adult                      | ARV120 (Green)  |
| EC-3830LK/P/FK2 Adult                 | ARV120 (Green)  |
| EC-380LKP/FKP/FK2P/MK2PAdult          | ARV120 (Green)  |
| EC-38LV Adult                         | ARV120 (Green)  |
| EC-3490FK Pediatric HD                | ARV130 (Purple) |
| EC-3490 FI Pediatric HD               | ARV130 (Purple) |
| EC-3490 MK/FK/LK/LI/TLI/TFI Pediatric | ARV130 (Purple) |
| EC-3470LK Pediatric                   | ARV130 (Purple) |
| EC-34-i10M/F/L Pediatric HD           | ARV130 (Purple) |
| EC-3490TLF/MI                         | Not compatible  |
| EC-2990LI                             | Not compatible  |

### STORTZ

|                     |                 |
|---------------------|-----------------|
| 13924-PKS Adult     | ARV120 (Green)  |
| 13925-PKS Adult     | ARV120 (Green)  |
| 13907-PKS Pediatric | ARV130 (Purple) |
| 13906-PKS Pediatric | ARV130 (Purple) |
| 13905-PKS Adult     | ARV140 (Orange) |
| 13904-PKS Adult     | ARV140 (Orange) |

\* The table provides an indication of ENDOCUFF VISION™ compatibility with colonoscopes. For confirmation of specific models, contact an Olympus representative.

Specifications, design, and accessories are subject to change without any notice or obligation on the part of the manufacturer. Olympus is a registered trademark of Olympus Corporation, Olympus America Inc., and/or their affiliates. ENDOCUFF VISION™ is a registered trademark of KeyMed (Medical & Industrial Equipment) Ltd.