

TC-C3
Compact Trolley



Compact Trolley

TC-C3

The TC-C3 Compact Trolley is designed as a compact trolley to be used in conjunction with Optera CV-170.

Compact

The compact design means the TC-C3 is ideal for use in the outpatient or doctor's office where space is limited.

Easy to Move

The Compact Trolley is easy to move and position due to its light weight and ergonomically designed handles at the front of the top tray.

Easy to Clean

The TC-C3 is easy to clean due to limited dirt traps and no exposed fittings.

Monitor Support

The TC-C3 is offered with a monitor mount that supports monitors up to 26" with unrestricted swivel and tilt movement, allowing good viewing positions of the procedure image.

Large Top Tray

The top tray provides a clear working space at a comfortable height. The increased range of accessories enables the trolley to be utilized for a range of procedures.

Integrated Separation Transformer

The TC-C3 can be supplied with a separating transformer that is placed under the base, this means increased shelf space for equipment loading.











Accessories

TC-C3

MAJ-2146

Sliding Keyboard Tray



MAJ-1686

Drawer Unit



MAJ-1656

Footswitch Holder



MAJ-2145 ENT Scope Hanger



MAJ-1684

Scope Pole Kit



MAJ-1685

IV Bag Hanger



MAJ-1639

CO₂ Cylinder Holder (140 mm)



MAJ-1658

Universal Storage Container



MAJ-2148

Rear Panel



MAJ-1683

LCD Monitor Mount



MAJ-1665

Camera Head Holder



MAJ-110

Jar Mounting Kit TC



Specifications	
Dimensions	1,110 mm (H) × 563.9 mm (D) × 489.5 mm (W)
	Height to top tray work surface - 1,100 mm
	Height with 26" LCD monitor in highest position – 1,945 mm
Weight	Maximum 82 kg (unladen) including fitted transformer

Voltage		
The TC-C3 can be supplied with a separating transformer based on one of the following voltage options:		
MAJ-2205	Transformer 110-120 V	
MAJ-2206	Transformer 220-240 V	

Trolley Construction		
Columns	Extruded profiled aluminum	
Shelves	Mild steel	
Top Tray	Injection-moulded ABS	

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

