

ENT Ordering Overview Rhinology

Selected Product Portfolio



CELON System

Procedure-Oriented, Safe and Efficient Solutions for ENT Specialists

CELON ELITE (ESG-200) Power Control Unit



- Is based on radio frequency and allows precise and safe tissue ablation for a variety of ENT indications.
- Can be used to perform standard electrosurgery in the ENT office or operating room.

Order Number	Product Description	Specifications
WA90001A	CELON ELITE, type E	Monopolar output socket for international/Erbe plugs (incl. double-pedal foot switch)
WA90002A	CELON ELITE, type B	Monopolar output socket for international/Bovie plugs (incl. double-pedal foot switch)
WA95621A	EU power cable	4.5 m with angled plug (many European countries, type E/F)
WA95622A	US power cable	4.5 m with angled plug (USA, Canada and other countries, type B)
WA95623A	UK power cable	4.5 m with angled plug (United Kingdom and other countries, type G)
WB50402W	Double-pedal foot switch	

Compatibility



Bipolar

For bipolar HF cables



International Standard

2-pin plug

► 4 mm diameter



Erbe Standard

1-pin coaxial plug

► 8/4 mm diameter



Monopolar

For monopolar HF cables

Please ensure that you connect the banana plug only to the rightmost receptacle.



International Standard

3-pin plug

► 4 mm diameter



Erbe Standard

1-pin coaxial plug

► 9/5 mm diameter



Banana Plug

1-pin plug

► 4 mm diameter



Neutral



Split and Nonsplit Neutral Electrodes

· 2-socket plugs

· 2.5 mm socket diameter

· 10 mm socket spacing

CELON System

Procedure-Oriented, Safe and Efficient Solutions for ENT Specialists



Bipolar RFITT Applicators



For minimally invasive, impedance-controlled tissue ablation when treating turbinate and tonsil hypertrophy, snoring and mild sleep apnea.

Order Number	Product Description	Connection	Unit
WB990004	① CelonProBreath	Celon LabENT	5 pcs./box
WB990095	② CelonProSleep Plus	Celon LabENT	5 pcs./box
WB990007	① CelonProBreath	Celon Elite	5 pcs./box
WB990008	② CelonProSleep Plus	Celon Elite	5 pcs./box

ProCut Accessories




- For monopolar cutting with the advantages of bipolar technology.
- Grasping forceps available in different shapes with various tip designs.

Order Number	Product Description	Shape	Tip Design	Length	Unit
WB990177	① Celon ProCut electrode				5 pcs./box
WB990278	② Celon Elite ProCut handpiece				1 pc./box
WB990203	③ Celon ProCut forceps	Angled	Grasping teeth	180 mm	1 pc./box
WB990204	④ Celon ProCut forceps	Bayonet, slim	Grasping teeth	180 mm	1 pc./box
WB990205	⑤ Celon ProCut forceps	Angled, slim	Grasping teeth	180 mm	1 pc./box
WB990201	⑥ Celon ProCut forceps	Straight	Atraumatic	145 mm	1 pc./box
WB990201	⑦ Celon ProCut forceps	Straight	Atraumatic	145 mm	1 pc./box

CELON System

Procedure-Oriented, Safe and Efficient Solutions for ENT Specialists

Bipolar Coagulation Forceps



- For effective vessel sealing.
- Available in different shapes.

Order Number	Product Description	Shape	Length	Tip Width	Tip Length	Tip Angle	Unit
WB990038	① Bipolar coagulation forceps	Angled	180 mm	2 mm	8,5 mm	45°	1 pc./box
WB990041	② Bipolar coagulation forceps	Bayonet	200 mm	1 mm	8,5 mm	0°	1 pc./box
WB990037	Bipolar coagulation cable, Valleylab plug, 3 m						1 pc./box


Other

Order Number	Product Description	Weight	Contact area	Unit
WB990200	Monopolar handpiece, hand and foot switch activation			
WA99512C	Neutral electrode, for newborns, with cable	< 5 kg	40 cm ²	15 pcs./box
WA99513A	Neutral electrode, for children, without cable	5-15 kg	70 cm ²	100 pcs./box
WA99513C	Neutral electrode, for children, with cable	5-15 kg	70 cm ²	15 pcs./box
WA99514A	Neutral electrode, for adults, without cable	> 15 kg	140 cm ²	100 pcs./box
WA99514C	Neutral electrode, for adults, with cable	> 15 kg	140 cm ²	15 pcs./box
WA99516A	Neutral electrode, universal size for children/adults without cable	> = 5 kg	110 cm ²	100 pcs./box
WA99516C	Neutral electrode, universal size for children/adults with cable	> = 5 kg	110 cm ²	15 pcs./box
WA99510A	Cable for neutral electrode, suitable for manual and automatic reprocessing and autoclaving, 4.5 m			
FSW991.067.003	Celon LabENT foot switch			

CELON System

Procedure-Oriented, Safe and Efficient Solutions for ENT Specialists

HF Accessories



- ① Needle electrode
- ② Blade electrode
- ③ Ball electrode
- ④ Loop electrode
- ⑤ Electrosurgical pencil
- ⑥ Suction coagulator

The Olympus electro-surgical portfolio offers a wide range of electro-surgical instruments, including disposable monopolar electrodes, pencils and suction coagulators, reusable bipolar forceps and a number of different accessories. Featuring a lightweight, ergonomic design and flexible cable, the Olympus electro-surgical pencils are preloaded with high quality and economical nonstick (PTFE) or stainless steel blade electrodes. To provide even more choice and value for all of the surgeon's needs, a variety of nonstick and stainless steel electrodes in different lengths and tip styles is additionally offered.

System usable with the ESG-200 and other standard high frequency generators.





CELON System

Procedure-Oriented, Safe and Efficient Solutions for ENT Specialists

HF Accessories

Ordering Information				
Order Number	Product Description	Size/Length/Diameter	Options	Pieces
Disposable Monopolar Electrosurgical Electrodes				
WA98510E	Electrode, blade	2.75"	Coated	10 pcs./box
WA98511E	Electrode, blade	2.75"	Non-coated	10 pcs./box
WA98520E	Electrode, blade	2.75"	Coated, extended insulation	10 pcs./box
WA98521E	Electrode, blade	2.75"	Non-coated, extended insulation	10 pcs./box
WA98530E	Electrode, blade	4"	Coated	10 pcs./box
WA98531E	Electrode, blade	4"	Non-coated	10 pcs./box
WA98540E	Electrode, blade	4"	Coated, extended insulation	10 pcs./box
WA98541E	Electrode, blade	4"	Non-coated, extended insulation	10 pcs./box
WA98550E	Electrode, blade	6"	Coated	10 pcs./box
WA98551E	Electrode, blade	6"	Non-coated	10 pcs./box
WA98560E	Electrode, blade	6"	Coated, extended insulation	10 pcs./box
WA98561E	Electrode, blade	6"	Non-coated, extended insulation	10 pcs./box
WA98610E	Electrode, needle	2.75"	Coated	10 pcs./box
WA98611E	Electrode, needle	2.75"	Non-coated	10 pcs./box
WA98620E	Electrode, needle	2.75"	Coated, extended insulation	10 pcs./box
WA98621E	Electrode, needle	2.75"	Non-coated, extended insulation	10 pcs./box
WA98630E	Electrode, needle	4"	Coated	10 pcs./box
WA98631E	Electrode, needle	4"	Non-coated	10 pcs./box
WA98640E	Electrode, needle	4"	Coated, extended insulation	10 pcs./box
WA98641E	Electrode, needle	4"	Non-coated, extended insulation	10 pcs./box
WA98650E	Electrode, needle	6"	Coated	10 pcs./box
WA98651E	Electrode, needle	6"	Non-coated	10 pcs./box
WA98660E	Electrode, needle	6"	Coated, extended insulation	10 pcs./box
WA98661E	Electrode, needle	6"	Non-coated, extended insulation	10 pcs./box
WA98830E	Electrode, ball	2.75", 5 mm	Coated	10 pcs./box
WA98840E	Electrode, ball	2.75", 5 mm	Non-coated	10 pcs./box
WA98810E	Electrode, ball	5.38", 5 mm	Coated	10 pcs./box
WA98820E	Electrode, ball	5.38", 5 mm	Non-coated	10 pcs./box
WA98950E	Electrode, loop	10 x 10 mm		10 pcs./box
WA98951E	Electrode, loop	15 x 10 mm		10 pcs./box
WA98952E	Electrode, loop	20 x 10 mm		10 pcs./box
WA98953E	Electrode, loop	20 x 15 mm		10 pcs./box
Disposable Monopolar Suction Coagulators				
WA98410C	Suction coagulator	8 Fr	Hand switch	15 pcs./box
WA98420C	Suction coagulator	10 Fr	Hand switch	15 pcs./box
WA98430C	Suction coagulator	12 Fr	Hand switch	15 pcs./box
WA98411C	Suction coagulator	8 Fr	Foot control	15 pcs./box
WA98421C	Suction coagulator	10 Fr	Foot control	15 pcs./box
WA98431C	Suction coagulator	12 Fr	Foot control	15 pcs./box



CELON System

Procedure-Oriented, Safe and Efficient Solutions for ENT Specialists

HF Accessories

Ordering Information				
Order Number	Product Description	Size/Length/Diameter	Options	Pieces
Disposable Monopolar Electrosurgical Pencils				
WA98210B	Electrosurgical pencil	3 m cable	Push button with coated blade	50 pcs./box
WA98211B	Electrosurgical pencil	3 m cable	Push button with non-coated blade	50 pcs./box
WA98220B	Electrosurgical pencil	3 m cable	Push button with coated blade, holster	30 pcs./box
WA98221B	Electrosurgical pencil	3 m cable	Push button with non-coated blade, holster	30 pcs./box
WA98230B	Electrosurgical pencil	3 m cable	Push button with coated high insulation blade, holster	30 pcs./box
WA98231B	Electrosurgical pencil	4.5 m cable	Push button with coated blade	50 pcs./box
WA98210R	Electrosurgical pencil	3 m cable	Rocker switch with coated blade	50 pcs./box
WA98211R	Electrosurgical pencil	3 m cable	Rocker switch with non-coated blade	50 pcs./box
WA98220R	Electrosurgical pencil	3 m cable	Rocker switch with coated blade, holster	30 pcs./box
WA98221R	Electrosurgical pencil	3 m cable	Rocker switch with non-coated blade, holster	30 pcs./box
WA98230R	Electrosurgical pencil	3 m cable	Rocker switch with coated high insulation blade, holster	30 pcs./box
WA98231R	Electrosurgical pencil	4.5 m cable	Rocker switch with coated blade	50 pcs./box
WA98210F	Electrosurgical pencil	3 m cable	Foot control with coated blade	50 pcs./box
WA98211F	Electrosurgical pencil	3 m cable	Foot control with non-coated blade	50 pcs./box
WA98220F	Electrosurgical pencil	3 m cable	Foot control with coated blade, holster	30 pcs./box
WA98221F	Electrosurgical pencil	3 m cable	Foot control with non-coated blade, holster	30 pcs./box
Disposable Smoke Evacuation				
WA98192S	OR-VAC monopolar electrosurgical pencil adapter	22 mm		20 pcs./box
Disposable Neutral Electrodes				
WA99512C	Neutral electrode, for newborns		With cable	15 pcs./box
WA99513A	Neutral electrode, for children		Without cable	100 pcs./box
WA99513C	Neutral electrode, for children		With cable	15 pcs./box
WA99514A	Neutral electrode, for adults		Without cable	100 pcs./box
WA99514C	Neutral electrode, for adults		With cable	15 pcs./box
WA99516A	Neutral electrode, universal size		Without cable	100 pcs./box
WA99516C	Neutral electrode, universal size		With cable	15 pcs./box
WA99510A	Cable, for neutral electrode, suitable for manual and automatic reprocessing and autoclaving*	4.5 m		1 pc./box
Reusable Bipolar Electrosurgical Devices				
WB990038	Bipolar forceps angled	180 mm		1 pc./box




DIEGO ELITE Multidebrider

Ordering Information – Blades and Burrs

Sinus Application Blades, Straight

Standard




For advanced resection of soft tissue.					
Order Number	Product Description	Diameter	Shape	Tip	Unit
EGSB4000SA	Shaver blade, type A	4 mm	Straight	Standard	5 pcs./box




For standard soft-tissue removal.					
Order Number	Product Description	Diameter	Shape	Tip	Unit
EGSB2900SC	Shaver blade	2.9 mm	Straight	Serrated	5 pcs./box
EGSB4000SS	Shaver blade	4 mm	Straight	Serrated	5 pcs./box

Bipolar



For resection of soft tissue and bipolar hemostasis.					
Order Number	Product Description	Diameter	Shape	Tip	Unit
EGBB4000SC	Bipolar blade	4 mm	Straight	Serrated	5 pcs./box

Monopolar



For resection of soft tissue and monopolar hemostasis.					
Order Number	Product Description	Diameter	Shape	Tip	Unit
EGMB4000SC	Monopolar blade	4 mm	Straight	Serrated	5 pcs./box

DIEGO ELITE Multidebrider

Ordering Information – Blades and Burrs

Sinus Application Blades, Curved

Standard



For standard soft-tissue removal in the different nasal cavities.

Order Number	Product Description	Diameter	Angulation	Tip	Unit
EGSB4015SC	Shaver blade	4 mm	15°	Serrated	5 pcs./box
EGSB2940SC	Shaver blade	2.9 mm	40°	Serrated	5 pcs./box
EGSB4040SC	Shaver blade	4 mm	40°	Serrated	5 pcs./box
EGSB4040RC	Shaver blade, rotatable window	4 mm	40°	Serrated	5 pcs./box
EGSB4060SC	Shaver blade	4 mm	60°	Serrated	5 pcs./box
EGSB4060RC	Shaver blade, rotatable window	4 mm	60°	Serrated	5 pcs./box
EGSB2975SC	Shaver blade	2.9 mm	75°	Serrated	5 pcs./box
EGSB4075RC	Shaver blade	2.9 mm	75°	Serrated	5 pcs./box
EGSB4090RC	Shaver blade, rotatable window	4 mm	90°	Serrated	5 pcs./box

Bipolar



For standard soft-tissue removal and bipolar hemostasis.

Order Number	Product Description	Diameter	Angulation	Tip	Unit
EGBB4040SC	Bipolar blade	4 mm	40°	Serrated	5 pcs./box

DIEGO ELITE Multidebrider

Ordering Information – Blades and Burrs

Turbinate Application Blades

Standard



2 mm turbinator blades for advanced submucosal resection.

Order Number	Product Description	Diameter	Angulation	Tip	Unit
EGSB2000SA	Turbinator shaver blade	2 mm	Straight	SMR	5 pcs./box


Bipolar



2 mm turbinator blades for advanced submucosal resection and bipolar hemostasis.

Order Number	Product Description	Diameter	Angulation	Tip	Unit
EGBB2000SA	Turbinator bipolar blade, type A	2 mm	Straight	SMR	Straight

Burrs



Selection of a wide range of different burrs for advanced FESS procedures.

Order Number	Product Description	Diameter	Angulation	Unit
EGBUR3015DB	① Diamond bullet burr	3 mm	15°	5 pcs./box
EGBUR3070DB	① Diamond bullet burr	3 mm	70°	5 pcs./box
EGBUR4015DT	② Diamond taper burr	4 mm	15°	5 pcs./box
EGBUR4070DT	② Diamond taper burr	4 mm	70°	5 pcs./box
EGBUR4515BB	③ Diamond ball burr	4.5 mm	15°	5 pcs./box
EGBUR4215FB	④ Fluted barrel burr	4.2 mm	15°	5 pcs./box
EGBUR4260FB	④ Fluted barrel burr	4.2 mm	60°	5 pcs./box

DIEGO ELITE Multidebrider

Ordering Information – Blades and Burrs

Tonsils and Adenoid Application Blades

Standard, Malleable



Malleable blades for tonsil and adenoid resection, bendable up to 40° (Convex angle).

Order Number	Product Description	Diameter	Angulation	Tip	Unit
EGMS4000SC	Malleable blade, T/A	4 mm	0-40°	Serrated	5 pcs./box


Bipolar, Curved



For tonsil and adenoid resection and bipolar hemostasis with convex 40-degree angulation.

Order Number	Product Description	Diameter	Angulation	Tip	Unit
EGBB4040XC	Bipolar blade, convex, T/A	4 mm	40°	Serrated	5 pcs./box

Monopolar, Malleable



Malleable blade for tonsil and adenoid resection with monopolar hemostasis, bendable up to 40° (Convex angle).

Order Number	Product Description	Diameter	Angulation	Tip	Unit
EGMM4000SC	Malleable monopolar blade, T/A	4 mm	0-40°	Serrated	5 pcs./box

Laryngeal Application Blades

Superficial



Laryngeal blade for removal of superficial laryngeal lesions.

Order Number	Product Description	Diameter	Angulation	Unit
EGLB2922SS	Laryngeal blade, superficial	2.9 mm	Straight	5 pcs./box
EGLB2927SS	Laryngeal blade, superficial	2.9 mm	Straight	5 pcs./box

Debulking



Laryngeal blade for tumor debulking.


Order Number	Product Description	Diameter	Angulation	Unit
EGLB4022SD	Laryngeal blade, debulking	4 mm	Straight	5 pcs./box
EGLB4027SD	Laryngeal blade, debulking	4 mm	Straight	5 pcs./box



DIEGO ELITE Multidebrider

Ordering Information – System and Accessories

System Compatibility Chart



DIEGO ELITE System

Order Number	Description
EGMDCONS100	Console
EGMDFS100	Foot switch
EGMDHP100A	MULTIDEBRIDER handpiece, angled
EGMDHP100S	MULTIDEBRIDER handpiece, straight
EGMDHPTRAY	MULTIDEBRIDER handpiece tray
EGTS100S	MULTIDEBRIDER tube set, standard
EGTS101DC	MULTIDEBRIDER tube set, declog

DIEGO ELITE Accessories

Order Number	Description	Unit
EGTS100S	MULTIDEBRIDER tube set, standard	5 pcs./box
EGTS101DC	MULTIDEBRIDER tube set, declog	5 pcs./box
K10021893	Suction module tube set, suction/filter	5 pcs./box
K10021894	Suction module tube set, suction/filter	5 pcs./box
K10026130	Suction module tube set, suction/filter	5 pcs./box
EGBEMTS101DC	Bracket, EM sensor	5 pcs./box
K10026145	KV-6 suction pump UK-ROW	1 pcs./box
K10026146	KV-6 suction pump EU	1 pcs./box
K10021893	Jars for liners 1.5liters	1 pcs./box
K10021894	Disposable liner 1.5 liters	40 pcs./box
K10026310	Suction tube 7mm x 12 mm connectors, 60cm (from Jars to generator)	1 pcs./box
K10037470	Patient connecting tube, 2m	15 pcs./box



InstaClear

Clear View of the Surgical Site, Anytime


InstaClear System



Lens cleaner system for instant cleaning of rigid endoscopes during surgeries without having to remove it from the nasal cavity.

Order Number	Product Description	Unit
EGLCCONS100	① InstaClear console	1 pc./box
EGLCFP100	② InstaClear foot switch	1 pc./box
EGLCTS100S	③ InstaClear tube set	5 pcs./box

Olympus InstaClear Sheaths



Compatible with Olympus 4 mm HD and ULTRA rigid endoscopes. Available with different working angles, as well as top or bottom light-cable positions.

Order Number	Product Description	Compatible Sinuscope	Light Position	Viewing Angle	Unit
EGLCS1500UNOL	Olympus InstaClear sheath	WA96200A	Universal	0°	5 pcs./box
EGLCS1530BTOL	Olympus InstaClear sheath	WA96201A	① Bottom	30°	5 pcs./box
EGLCS1530TPOL	Olympus InstaClear sheath	WA96202A	② Top	30°	5 pcs./box
EGLCS1545BTOL	Olympus InstaClear sheath	WA96203A	① Bottom	45°	5 pcs./box
EGLCS1545TPOL	Olympus InstaClear sheath	WA96204A	② Top	45°	5 pcs./box
EGLCS1570BTOL	Olympus InstaClear sheath	WA96205A	① Bottom	70°	5 pcs./box
EGLCS1570TPOL	Olympus InstaClear sheath	WA96206A	② Top	70°	5 pcs./box

Olympus ULTRA InstaClear Sheaths


Order Number	Product Description	Compatible Telescope	Light Position	Viewing Angle	Unit
EGLCS4K00UNOL	Olympus ULTRA InstaClear sheath	WA4KS400	Universal	0°	5 pcs./box
EGLCS4K30BTOL	Olympus ULTRA InstaClear sheath	WA4KS430	① Bottom	30°	5 pcs./box
EGLCS4K30TPOL	Olympus ULTRA InstaClear sheath	WA4KS431	② Top	30°	5 pcs./box
EGLCS4K45BTOL	Olympus ULTRA InstaClear sheath	WA4KS445	① Bottom	45°	5 pcs./box
EGLCS4K45TPOL	Olympus ULTRA InstaClear sheath	WA4KS446	② Top	45°	5 pcs./box
EGLCS4K70BTOL	Olympus ULTRA InstaClear sheath	WA4KS470	① Bottom	70°	5 pcs./box
EGLCS4K70TPOL	Olympus ULTRA InstaClear sheath	WA4KS471	② Top	70°	5 pcs./box



InstaClear

Clear View of the Surgical Site, Anytime

Storz InstaClear Sheaths




Compatible with 4 mm Storz HOPKINS sinuscope. Available with different working angles, as well as top or bottom light-cable positions.

Order Number	Product Description	Compatible HOPKINS Sinuscope	Light Position	Viewing Angle	Unit
EGLCS1800UNST	Storz InstaClear sheath	7230AA	Universal	0°	5 pcs./box
EGLCS1830BTST	Storz InstaClear sheath	7230BA	① Bottom	30°	5 pcs./box
EGLCS1830TPST	Storz InstaClear sheath	7230BVA	② Top	30°	5 pcs./box
EGLCS1845BTST	Storz InstaClear sheath	7230FA	① Bottom	45°	5 pcs./box
EGLCS1845TPST	Storz InstaClear sheath	7230FVA	② Top	45°	5 pcs./box
EGLCS1870BTST	Storz InstaClear sheath	7230CA	① Bottom	70°	5 pcs./box
EGLCS1870TPST	Storz InstaClear sheath	7230CVA	② Top	70°	5 pcs./box

Sinuscopes

A Powerful Solution for Sinus Surgery


HD Sinuscopes



HD sinuscopes with outstanding image quality. Available with different working angles, as well as top or bottom light-cable positions. Autoclavable, 146 mm length.

Order Number	Product Description	Diameter	Light position	Viewing Angle	Unit
WA96200A	HD sinuscope	4 mm	Universal	0°	1 pc./box
WA96201A	HD sinuscope	4 mm	① Bottom	30°	1 pc./box
WA96202A	HD sinuscope	4 mm	② Top	30°	1 pc./box
WA96203A	HD sinuscope	4 mm	① Bottom	45°	1 pc./box
WA96204A	HD sinuscope	4 mm	② Top	45°	1 pc./box
WA96205A	HD sinuscope	4 mm	① Bottom	70°	1 pc./box
WA96206A	HD sinuscope	4 mm	② Top	70°	1 pc./box

ULTRA Telescopes



The ULTRA telescopes are optimized for usage with high-resolution camera systems (HD and 4K) to provide state-of-the-art image quality to the customer. Available with different working angles, as well as top or bottom light-cable positions. Autoclavable, 176.5 mm length.

Order Number	Product Description	Diameter	Light position	Viewing Angle	Unit
WA4KS400	ULTRA telescope	4 mm	Universal	0°	1 pc./box
WA4KS430	ULTRA telescope	4 mm	① Bottom	30°	1 pc./box
WA4KS431	ULTRA telescope	4 mm	② Top	30°	1 pc./box
WA4KS445	ULTRA telescope	4 mm	① Bottom	45°	1 pc./box
WA4KS446	ULTRA telescope	4 mm	② Top	45°	1 pc./box
WA4KS470	ULTRA telescope	4 mm	① Bottom	70°	1 pc./box
WA4KS471	ULTRA telescope	4 mm	② Top	70°	1 pc./box

Sinuscope Accessories

Order Number	Product Description	Size/Diameter	Unit
WA05990A	Wire mesh tray, for 2 telescopes	446 × 49 × 88 mm	1 pc./box
WA03300A	Light-guide cable, 3 m	2.8 mm	1 pc./box
WA03310A	Light-guide cable, 3 m	4.25 mm	1 pc./box

CELERIS

The New Portable Sinus Debrider

CELERIS System



- Cost effective and portable option for basic nasal procedures such as polypectomies, mini FESS and turbinate reduction.
- Comprehensive portfolio with four options to cover a wide range of surgical needs
- One device replaces multiple blades
- Manages incidental bleeding without instrument exchange when connected to a bipolar electro surgical generator
- Quick and easy plug and perform setup

Ordering Information

Order Number	Product Description	Specifications
CELERIS portfolio		
EGDSD4000SA	Celeris Disposable Debrider. 4mm standard	4mm straight standard blade
EGDSD40MMSC	Celeris Disposable Debrider. 4mm malleable	4mm straight standard blade, Malleable up to 60°
EGDSD4000BA	Celeris Disposable Debrider. 4mm bipolar	4mm straight blade with Bipolar energy (must be connected to a compatible bipolar electro surgical generator to activate bipolar energy)
EGDSD4000BA	Celeris Disposable Debrider. 2mm bipolar	2mm straight blade with Bipolar energy for Turbinate SMR (must be connected to a compatible bipolar electro surgical generator to activate bipolar energy)
EGDSDPP100	Celeris Reusable Power Pack	Power Pack for operating and controlling the disposable debrider
CELERIS optional accessories		
EGLCFP100	Foot pedal	Foot switch for foot activation
WA90001A	CELON ELITE ESG 200 & double pedal foot switch	Radio frequency based power control unit for RFITT and standard electro surgical procedures. Controls CELERIS bipolar blades
K10026145	KV-6 Suction Pump UK-ROW	Suction system for endoscopic aspiration and for non endoscopy applications
K10026146	KV-6 suction pump EU	
K10026128	Reusable Jars 0.25 litres (1pcs/box)	
K10026127	Reusable lids for 0.25 reusable jars (1pcs/box)	
K10026310	Suction tube 7mm x 12 mm connectors, 60cm (from Jars to generator) (1pcs)	
K10037470	Patient connecting tube, 2m (15pcs/box)	

VISERA ELITE III

The New Endoscopic Visualization Platform for ENT

VISERA ELITE III System

Visera Elite III allows you to discover 4K and 3D imaging and well proven Narrow Band Imaging (NBI) all supported in one system.

- Advanced imaging for improved surgical outcome
- Narrow Band Imaging (NBI) for improved diagnostic outcome
- Optimized picture with Flexible Rhinolaryngoscopes
- Enables 4K rigid endoscopy and full HD flexible endoscopy



Ordering Information

Order Number	Product Description	Specifications/features
VISERA ELITE III portfolio		
E0497746	LMD-XH320ST 32" LCD Monitor	4K Resolution, 2D & 3D Compatibility, A.I.M.E (Advanced Image Multiple Enhancer), 12G SDI, input
E0497747	LMD-XH550ST 55" LCD Monitor	
N608160	OEV321UH 4K Monitor	Enables full 4k image quality and excellent color reproduction
E0421076	LMD-X27MD 27" LCD Sony Monitor	4K Resolution
N6158050	OTV-S700 Video Processor	4K and NBI, 3D can easily be added via software upgrades, HDR, 12G SDI output
N6157950	CLL-S700 LED Light Source	Automatic brightness control, HDR support, LED lamp
N6154750	CH-S700 XZ EA 4K Camera Head	Autoclave compatibility, Continuous Auto Focus (CAF) & Extended Depth of Field (EDOF) Ergonomic
N6189550	IR Software Upgrade Activation USB Memory Key	Simply update as upgrade and innovation become available



VISERA ELITE II

A Unique World of Observation Possibilities

VISERA ELITE II System



Visera Elite II is an all in one video system optimized for ENT examinations and office based procedures. It is compatible with a wide range of rigid and flexible endoscopes, camera heads.

- Excellent natural color reproduction with enhanced imaging
- Intuitive handling and easy set up
- NBI and two modes of IR observation
- Compact one box system to save costs, time and space

Ordering Information		
Order Number	Product Description	Specifications
OTV S200 Video System Center		
N5426850	OTV-S200 Video System	2D processor and light source
N5766550	Full HD Camera Head Autoclavable	Excellent image quality and color reproduction, autoclavable
E0421076	Sony LMD-X270MD 27" 4K LCD Monitor	Enables full HD upscaling for outstanding image quality
E0421077	NDS Radiance Ultra 27" HD Monitor	Full HD monitor



nCare

Secure and Intuitive Design to Support Health Care Teams

nCare System



Uncover the remarkable benefits of a tailored centralized documentation system thoughtfully designed for ENT practices Effortlessly guide patients throughout their entire care journey, from initial evaluation to post treatment follow up, nurturing uninterrupted continuity of care, optimal treatment outcomes, and enhanced patient satisfaction.

nCare medical recorder can be further integrated into a centralized solution for documentation and archiving of medical information and medical collaboration* This technology integrates seamlessly with existing workflows, connecting health care teams and enabling them to focus on delivering outstanding patient care now and with changing demands in the future.

Ordering Information

Order Number	Product Description
nCare portfolio	
U9000407	nCare medical recorder, Single channel, SDI
U9000409	nCare medical recorder, Dual channel, 2xSDI
nCare accessories	
U9000698	ELO Commercial Grade Touchscreen, 15 inch
E0420542	Canvys Medical Grade Touchscreen, 15.6 inch
K10030210	Mounting Arm for touchscreen (MAJ-2167)



*For further information please contact your local sales representative.