OLYMPUS

KV-6 Medical Virus Filters





Medical Virus Filters

Be on the safe side and reduce the risk of viral spread

Electrical, mobile suction pumps, equipped with filters and accessories help to prevent viral spread in hospitals and care centres. We offer one complete infection control filter portfolio and fulfil new and strictest requirements (ISO 10079).

Safe

- \cdot Lab-tested high virus retention grades of **99.99984** % for particles larger than 0.03 μm
- · Filters all types of Coronavirus: 2019-nCoV, SARS-CoV, HCoV NL 63 and HKUI
- · Also filters other virus and bacteria types

High Performance

- · Dedicated, high-performance medical membranes
- · Low air resistance due to the large surface area


Flexible

- · Olympus' disposable virus filters are seamlessly compatible with Olympus KV-6 for endoscopic, airway and surgical suction
- · Custom-configure your suction system to minimize the risk of cross-contaminations



Medical Virus Filters

Filtration Efficiency - µm (Micrometer)





Filter Replacement

Filter should be replaced -

- in the event of overflow or clogging, or
- weekly with single patient use, or
- daily with multiple patient use

To detect a clogged filter, perform a filter test before each patient use, always complying with in-house directives on hygiene.

Medical Virus Filters

Ordering Information Medical Virus Filters		
K10027321	Disposable bacteria / odour filter	10 pack
K10027322	Disposable bacteria filter	10 pack
K10036997	Disposable viral / odour filter	10 pack
K10036992	Disposable viral filter	10 pack

As medical knowledge is constantly growing, technical modifications or changes of the product design, product specifications and accessories may occur.



For a complete listing of sales and distribution locations visit: www.olympus.com