

Irrigation Tubing for use with with CO₂ and Air



Irrigation Tubing with CO₂ / Air

Olympus irrigation tubing for use with CO₂ or Air are the 3-in-1 solution for insufflation, lens cleaning and irrigation from one water source.

Improve your cost per procedure by combining your existing insufflation / lens flush and irrigation tubes

- Improved irrigation flow vs existing systems
- 24 hour use
- Innovative bottle cap offers maximum seal with improved handling
- Dual connector design prevents leakage over multiple procedural changes
- Decreased set-up time - all in one package
- To be used with commercially available disposable water bottles
- Product requires MAJ-1652 Auxiliary Channel Adaptor
- CO₂ version requires connection tube for UCR: Single-Day-Use Contu CO₂ connection tube or reusable gas tube (MAJ-1741, MAJ-1742 & MAJ-1816)

Two piece bottle cap

No more twisted tubes



Ergonomic dual connector

Improved sealing and handling




Pump head markers

Easier placement





Traceability label

Prevent accidental re-use





MAJ-2207/MAJ-2208







Ordering Information

Irrigation Tubing with CO₂

Article Number	Article Name	Description	Note	Compatible Water Bottles (1000 ml)
 K10035002	MAJ-2207	Irrigation tubing with CO ₂ Cap A	Supplied sterile/ Single-Day-Use	Baxter UKF7114 / AKF7114 Baxter 2F7114 B. Braun R5000-01 Hospira 07139-09
 K10035003	MAJ-2208	Irrigation tubing with CO ₂ Cap B	Supplied sterile/ Single-Day-Use	Baxter UKF7114 / AKF7114 Fresenius Kabi 3158589 Fresenius Kabi Ampuwa 16-01-FAM10 Laboratoire Aguettant Versol 600500



CO₂ Tubing

Article Number	Article Name	Description	Note	Quantity - Unit/Box
 E0420575	Contu CO ₂	CO ₂ Connection Tube including removable UCR Adapter*	Single-Day-Use	10
 N2654500	MAJ-1741	UCR High Flow Gas Tube	Reusable	1
 N2654600	MAJ-1742	UCR Medium Flow Gas Tube	Reusable	1
 N3493330	MAJ-1816	UCR Low Flow Gas Tube	Reusable	1

* Distributed by: Olympus Europa SE & Co. KG
Manufactured by: MDD Company GmbH, Renglochstrasse 38, 6012 Obernau, Switzerland

Irrigation Tubing

MAJ-2209/MAJ-2210



Ordering Information

Irrigation Tubing with Air

Article Number	Article Name	Description	Note	Compatible Water Bottles (1000 ml)
82730 K10035004	MAJ-2209	Irrigation tubing with Air Cap A	Supplied sterile/ Single-Day-Use	Baxter UKF7114 / AKF7114 Baxter 2F7114 B. Braun R5000-01 Hospira 07139-09
82724 K10035005	MAJ-2210	Irrigation tubing with Air Cap B	Supplied sterile/ Single-Day-Use	Baxter UKF7114 / AKF7114 Fresenius Kabi 3158589 Fresenius Kabi Ampuwa 16-01-FAM10 Laboratoire Aguettant Versol 600500

As medical knowledge is constantly growing, technical modifications or changes of the product design, product specifications and accessories may be required.

OLYMPUS

OLYMPUS EUROPA SE & CO. KG
Postbox 10 49 08, 20034 Hamburg, Germany
Wendenstrasse 14-18, 20097 Hamburg, Germany
Phone: +49 40 23773-0, Fax: +49 40 233765
www.olympus-europa.com

For a complete listing of sales and distribution locations visit:
www.olympus.com