

# OLYMPUS

Your Vision, Our Future

## PVA SPONGE PRODUCTS

**Bacteriostatic PVA sponges**



## WINNER OF THE MEDICAL DESIGN EXCELLENCE AWARD IN 2000

Blu and Blu Glide sponge products take PVA packing to a new level. Available in all popular sizes for nasal, FESS, epistaxis and even otologic applications, Blu and Blu Glide provide benefits beyond what standard PVA sponge packs can offer. Both products contain a safe, patented, bacteriostatic compound, which inhibits the growth of infection and odour-causing microorganisms on the sponge.



## Bacteriostatic PVA sponges


In addition, Blu Glide possesses “glide surface” technology which helps prevent sticking upon removal of the sponge. This breakthrough patent-pending technology eliminates the problems associated with current laminated products.

Blu Glide “glide surface” is double-compressed for easier insertion and placement, and its microporous surface assures proper expansion with none of the cracking or separation associated with laminated layers.


If the answer is sponge packing, the choice is Blu or Blu Glide for an added measure of safety and patient comfort.

### Blu and Blu Glide are:

- Biocompatible
- Bacteriostatic
- Hydrophilic
- Micropore surface
- Durable yet impressively soft



*Leading competitor's PVA sponge (20x mag.)*



*Blu Glide PVA sponge (20x mag.)\**

### Blu and Blu Glide have been proven effective for inhibiting growth of the following troublesome bacteria and viruses on the foam for up to seven days:

- Staphylococcus aureus, including MRSA
- Staphylococcus epidermidis
- Bacillus subtilis
- Yersinia enterocolitica
- Serratia marcescens
- Vancomycin-resistant Enterococcus faecalis (VRE)
- Escherichia coli
- Vesicular stomatitis virus (VSV)
- Encephalomyocarditis virus (EMCV)

*\*Microscopic comparison of Blu Glide PVA foam to the leading product demonstrates smaller pore size and consistency of Blu and Blu Glide products.*

# PACKS

Blu and Blu Glide sponge products are bacteriostatic and biocompatible. Performance properties – such as expansion rate and softness – have been fine-tuned to assure ease of use and comfort. All Olympus ENT Blu and Blu Glide PVA nasal and sinus packs include hand-tied drawstrings for safety.

All dimensions reflect size in expanded (hydrated) state.

## Nasal packs

### Nasal pack, small bacteriostatic\*

Catalogue No.	Description	Length	Thickness	Height
EG70890830	Blu, bacteriostatic	4.5 cm	1.5 cm	2.0 cm
EG70890831	Blu, bacteriostatic, with glide technology	4.5 cm	1.5 cm	2.0 cm




### Nasal pack, large\*

Catalogue No.	Description	Length	Thickness	Height
EG70890832	Blu, bacteriostatic	8.0 cm	1.5 cm	2.0 cm
EG70890833	Blu, bacteriostatic, with glide technology	8.0 cm	1.5 cm	2.0 cm



### Nasal pack, small, with airway\*

Catalogue No.	Description	Length	Thickness	Height
EG70890834	Blu, bacteriostatic	4.5 cm	1.5 cm	2.0 cm
EG70890835	Blu, bacteriostatic, with glide technology	4.5 cm	1.5 cm	2.0 cm



### Nasal pack, large, with airway\*

Catalogue No.	Description	Length	Thickness	Height
EG70890836	Blu, bacteriostatic	8.0 cm	1.5 cm	2.0 cm
EG70890837	Blu, bacteriostatic, with Glide technology	8.0 cm	1.5 cm	2.0 cm



\* All nasal packs are packaged 1 pack per pouch, 10 pouches per box.

### D-style anatomical nasal pack\*

Catalogue No.	Description	Length	Thickness	Height
EG70890838	Blu, bacteriostatic	8.0 cm	1.5 cm	3.0 cm
EG70890839	Blu, bacteriostatic, with glide technology	8.0 cm	1.5 cm	3.0 cm




### D-style anatomical nasal pack, with airway\*

Catalogue No.	Description	Length	Thickness	Height
EG70890840	Blu, bacteriostatic	8.0 cm	1.5 cm	3.0 cm
EG70890841	Blu, bacteriostatic, with glide technology	8.0 cm	1.5 cm	3.0 cm



### D-style slim anatomical nasal pack\*

Catalogue No.	Description	Length	Thickness	Height
EG70890842	Blu, bacteriostatic	8.0 cm	0.75 cm	3.0 cm
EG70890843	Blu, bacteriostatic, with glide technology	8.0 cm	0.75 cm	3.0 cm



## Epistaxis packs

### Epistaxis pack, large\*

Catalogue No.	Description	Length	Thickness	Height
EG70890856	Blu, bacteriostatic	10.0 cm	1.5 cm	2.5 cm



### Epistaxis pack, small\*

Catalogue No.	Description	Length	Thickness	Height
EG70890857	Blu, bacteriostatic	8.5 cm	1.5 cm	2.5 cm




\* All nasal packs are packaged 1 pack per pouch, 10 pouches per box.

# PACKS

## Sinus packs

### K-style sinus pack\*

Catalogue No.	Description	Length	Thickness	Height
EG70890844	Blu, bacteriostatic	3.5 cm	0.9 cm	1.2 cm
EG70890845	Blu, bacteriostatic, with glide technology	3.5 cm	0.9 cm	1.2 cm




### K-style slim sinus pack\*

Catalogue No.	Description	Length	Thickness	Height
EG70890850	Blu, bacteriostatic	3.5 cm	0.6 cm	1.2 cm
EG70890851	Blu, bacteriostatic, with glide technology	3.5 cm	0.6 cm	1.2 cm




### B-style sinus pack\*

Catalogue No.	Description	Length	Thickness	Height
EG70890846	Blu, bacteriostatic	2.5 cm	0.9 cm	2.0 cm
EG70890847	Blu, bacteriostatic, with glide technology	2.5 cm	0.9 cm	2.0 cm




### B-style sinus pack\*

Catalogue No.	Description	Length	Thickness	Height
EG70890836	Blu, bacteriostatic	2.5 cm	0.6 cm	2.0 cm
EG70890837	Blu, bacteriostatic, with Glide technology	2.5 cm	0.6 cm	2.0 cm



### T-style sinus pack\*

Catalogue No.	Description	Length	Thickness	Height
EG70890836	Blu, bacteriostatic	4.5 cm	0.9 cm	2.0 cm
EG70890837	Blu, bacteriostatic, with Glide technology	4.5 cm	0.9 cm	2.0 cm



\* All sinus packs are packaged 2 packs per pouch, 10 pouches per box.

### T-style sinus pack\*

Catalogue No.	Description	Length	Thickness	Height
EG70890854	Blu, bacteriostatic	4.5 cm	0.6 cm	2.0 cm
EG70890855	Blu, bacteriostatic, with glide technology	4.5 cm	0.6 cm	2.0 cm



### Sinus kit

#### Sinus kit, bacteriostatic (for use following turbinectomy)\*

Catalogue No.	Description	Length	Thickness	Height
EG70890861	Blu, bacteriostatic	-	-	-
EG70890855	Blu, bacteriostatic, with glide technology	-	-	-



### Otology products

#### Blu ear wick

Catalogue No.	Description	Length	Thickness	Height
EG70890858	Bacteriostatic, 1 wick per pouch, 50 pouches per box	0.9 cm	1.5 cm	-
EG70890859	Bacteriostatic, 1 wick per pouch, 10 pouches per box	0.9 cm	2.4 cm	-



#### Blu ear packing with drawstring

Catalogue No.	Description	Length	Thickness	Height
EG70890860	Bacteriostatic, 1 wick per pouch, 10 pouches per box	1.2 cm	1.5 cm	-



\* All sinus packs are packaged 2 packs per pouch, 10 pouches per box.

# PVA SPONGE PRODUCTS

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

**OLYMPUS**

**OLYMPUS EUROPA HOLDING GMBH**

Postbox 10 49 08, 20034 Hamburg, Germany  
Wendenstrasse 14-18, 20097 Hamburg, Germany  
Phone: +49 (0)40 237 730, Fax: +49 (0)40 230 761  
[www.olympus-europa.com](http://www.olympus-europa.com)