

SafeSnap Polyp Traps

Single- and Multi-Chamber Polyp Traps to Meet Your Specimen Retrieval Needs




SafeSnap Polyp Traps

Designed for Efficient Sample Acquisition, Safe Specimen Handling and Convenient Installation



Indexing System
Allows correct localization and identification of multiple specimens.



Safety Screen
For filtering the specimen, helping to minimize exposure to bodily fluids.

Safety Screen
Designed to avoid loss of specimen and capture additional polyps.

Conical Adapter
Fits most suction tubes.



Safety
Individual chambers for sample separation and safe transportation, avoiding contamination.



Quality
Multiple chambers limit the need for specimen transfer during procedures, avoiding the drying out of the sample.



Time
Multiple chambers and safety screen facilitate the uninterrupted suction of various polyps, providing a smooth procedure workflow.



Fit
Conical adapter compatible with most suction tubes — one less worry!

SafeSnap Polyp Traps

SafeSnap 4+

Single-Use Polyp Trap with 4+ Chambers



SafeSnap 1

Single-Use 1 Chamber Polyp Trap



Ordering Information

Article Number	Article Name	Quantity	Description
E0420516	PT-01-04	25	SafeSnap 4+ Polyp Trap
E0420515	PT-01-01	25	SafeSnap 1 Polyp Trap

 www.et-catalogue.olympus.eu/



Olympus reserves the right of errors, modification and changes of the service and/or product offerings.

OLYMPUS

Distributed by: **OLYMPUS EUROPA SE & CO. KG**
Postbox 10 49 08, 20034 Hamburg, Germany
Wendenstrasse 14-18, 20097 Hamburg, Germany
Phone: +49 40 23773-0, Fax: +49 40 233765
www.olympus-europa.com

Manufactured by: **MDD Company GmbH**
Rengglochstrasse 38
6012 Obernau
Switzerland