

Single Use Air/Water, Suction & Biopsy Valve Kit

## Clever Shield

Dedicated Single-use Solution Streamlines Endoscopy Procedures



## Dedicated Single-use Solution Streamlines Endoscopy Procedures

### Performance Required for Endoscopy

The Clever Shield is a single-use valve kit from Olympus, the industry-leading endoscope manufacturer. Including air/water valve, suction valve and biopsy valve, the Clever Shield satisfies the performance at the Olympus-defined standard level as with the reusable valves, and provides secure air/water flow and suction rate with assured reliability and trusted security.



#### Safety and User-friendly

- Presterilized single-use valves
- Traceability assured by information stickers (sterilization information and lot no. )
- Allows immediate preparation by one-step-opening package
- Color usage differentiates from reusable valves



#### Single-use Air/Water Valve, Single-use Suction Valve

- Air/water flow and suction rate satisfy Olympus standards
- Color coded for identification just like reusable valves
- Silicone-free design
- Slightly rounded valve top supports air manipulation



#### Single-use Biopsy Valve

- Easy detachment by pressing down the lever
- Wide slit and unique material composition expected to contribute to smoother device manipulation
- Dual-chamber structure and small inner hole prevents leakage
- Also sold individually as MAJ-1555

Specifications			
Clever Shield/MAJ-2287		MAJ-1555	
Number of Units in Box	25	Number of Units in Box	20
Components	Single-use air/water valve Single-use suction valve Single-use biopsy valve	Components	Single-use biopsy valve
SINGLE USE AIR/WATER, SUCTION, & BIOPSY VALVE KIT (Sterile)		SINGLE USE BIOPSY VALVE (Sterile)	

As medical knowledge is constantly growing, technical modifications or changes of the product design, product specifications, accessories and service offerings may be required.