

**ENDO EYE 3D 10 MM**

**Video Telescope**



# ENDO-EYE 3D 10 MM

## 3D-Image Rotation without Loss of Horizon

- Change the view direction while maintaining horizontal orientation of images, enabling a continuous critical view in 3D

## Chip-on-the-Tip Technology

- Bright, clear, and natural 3D depth perception
- Focus-free handling; no manual focusing required

## Autoclave Compatibility

- Reduced costs compared with other sterilization methods

## Full Control During Surgery

- Three programmable remote control buttons

## Plug and Play

- Simple to set up and easy to use

### Ordering Information

Order Number	Outer Diameter	Direction of View	Field of View	Working Length	Max. Insertion Portion Width	Cable Length	Weight
WA50080A	10 mm	0°	67°	330 mm	10.07 mm	2780 mm	290g
WA50082A	10 mm	30°	67°	335 mm	10.20 mm	2780 mm	290g

### Specifications

<b>Reprocessing Methods</b>	Autoclavable, ETO, STERRAD, washer-disinfector		
<b>Ambient Conditions</b>	Operating conditions	Ambient temperature	10 to 35°C (50 to 95°F)
		Relative humidity	30 to 85%
	Storage conditions	Temperature	10 to 40°C (50 to 104°F)
		Relative humidity	30 to 75%
	Transport conditions	Temperature	-47 to 70°C (-53 to 158°F)
		Relative humidity	10 to 95%
<b>Classification of ME Equipment and ME Systems</b>	Protection against electric shock: Classification of applied parts	Type BF	
	Protection against harmful ingress of water or particulate matter	IPX7	

Specifications, design, and accessories are subject to change without any notice or obligation on the part of the manufacturer.

**OLYMPUS**

**OLYMPUS EUROPA SE & CO. KG**  
Postbox 10 49 08, 20034 Hamburg, Germany  
Wendenstrasse 14-18, 20097 Hamburg, Germany  
Phone: +49 40 23773-0, Fax: +49 40 233765  
www.olympus-europa.com

