


**MAJ-2294**

**Lead Antenna Unit**





## Features and Benefits

### Optimized Fit

The lead-type antenna fits various body shapes and is more suitable for bedridden patients.

### Easy and Reliable Reprocessing

The improved design eases reprocessing by making cable clearances accessible to mild detergent.\*



## Product Specifications

### Lead Antenna Unit (MAJ-2294)

<b>Weight</b>	200 g
<b>Dimensions</b>	87 (W) × 51 (H) × 15 (D) mm (without projection parts)
<b>Number of Leads</b>	8 antennae

### Antenna Lead Cover (MAJ-1470)

<b>Dimensions</b>	92 (W) × 66 (H) mm
<b>Quantity</b>	50 sheets/box

\* Please do not moisten the connection unit.

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.