



Single-use Biopsy Forceps

# EndoJaw

**Comprehensive Lineup Delivering Optimal Performance** 



## EndoJaw

### **Single-use Biopsy Forceps**

### **Deliver Optimal Performance**

Olympus offers an extensive range of single-use biopsy forceps. Featuring superior construction and delivering reliable performance, EndoJaw devices are designed to work seamlessly with a broad array of Olympus endoscopes.

#### **Comprehensive Lineup**

EndoJaw biopsy forceps are available in a wide variety of configurations, including alligator jaw and needle designs as well as a choice of diameters: 2.0 mm, 2.8 mm, 3.2 mm, and 3.7 mm.

# **Simplified Targeting** Proprietary swing jaw design mechanism facilitates **Alligator Jaw** the positioning of the forceps cups at biopsy sites with needle and allows for tangential approaches. 3.7 mm **Alligator Jaw** 3.2 mm **EndoJaw Large** Alligator jaw with needle 2.8 mm **Oval Jaw** with needle 2.8 mm **Oval Jaw** 2.0 mm



# EndoJaw

### **Specifications**

Specifications for Upper GI								
Article No.	Min. Working channel Ø	Cup type	Feature	Quantity	Working length	Sheath color		
FB-211K	2.0 mm _	Alligator Jaw -	-	- - - 20/box -	1550 mm	Grey		
FB-221K			With needle					
FB-231K		Oval Jaw -	-					
FB-241K			With needle					
FB-210K	2.8 mm	Alligator Jaw -	-					
FB-220K			With needle					
FB-230K		Oval Jaw -	-					
FB-240K			With needle					

pecifications	for Lower GI					
Article No.	Min. Working channel Ø	Cup type	Feature	Quantity	Working length	Sheath color
FB-210U		Alligator Jaw -	-		2300 mm	Blue
FB-220U	 2.8 mm		With needle			
FB-230U	Standard Capacity	Oval Jaw -	-			
FB-240U			With needle			
FB-215U	2.8 mm Large Capacity	Alligator la	-	- 20/DUX -		
FB-225U		Alligator Jaw -	With needle			
FB-235U		Oval Jaw -	-			
FB-245U			With needle			
FB-214U	3.2 mm -	Alligator Jaw	_	20/box		
FB-224U			With needle			
FB-234U		Oval Jaw	_			
FB-244U			With needle			
FB-212U	3.7 mm _	Alligator Jaw		5/box		
FB-222U			With needle			
FB-232U		Oval Jaw -	-			
FB-242U			With needle			

Specification for Respiratory								
Article No.	Min. Working channel Ø	Cup type	Feature	Quantity	Working length	Sheath color		
FB-211D	 2.0 mm - _	Alligator Jaw	-	- - 20/box -	1150 mm	Croon		
FB-221D			With needle					
FB-231D		Oval Jaw	-		1150 11111	Green		
FB-241D			With needle					

Specifications for Medical Thoracoscopy								
Article No.	Min. Working channel Ø	Cup type	Feature	Quantity	Working length	Sheath color		
FB-420K	2.8 mm	Alligator Jaw	With needle	20/box	1517 mm	Grey		

As medical knowledge is constantly growing, technical modifications or changes of the product design, product specifications, accessories and service offerings may be required.

