

ULTRASONIC TECHNOLOGY

Rapid cutting and reliable vessel sealing from a single multi-functional instrument.



Main features

- Wiper jaw design distributes pressure evenly along the jaw length for consistent sealing.
- Excellent grasper due to high grasping forces at the instrument tip.
- Excellent dissector due to fine tip design and a wide opening aperture.
- Rapid and reliable tissue cutting.
- Two handle types (inline and pistol grip) and four working lengths (10, 20, 35 and 45 cm) to meet a variety of surgical needs.
- Lower mist generation for consistently clearer surgical view.
- VAR and MAX handswitches on handpiece for surgeon control.

Generators and accessories

SONICBEAT

N3808660	USG-400	Ultrasonic generator
N3635730	TC-E400	Energy cart
N3808860	TD-SB400	SONICBEAT transducer
N3809130	MAJ-1869	Footswitch for SONICBEAT
N3809330	MAJ-1871	Communication cable 0.25 m
N3809430	MAJ-1872*	Communication cable 10 m
N3809530	MAJ-1873*	Adapter for UHI-2/3
N3809630	MAJ-1876	Docking fixture

*Required for the automatic mist and smoke evacuation function

SONICBEAT instruments

SONICBEAT

N3811530	SB-0545PC	SONICBEAT 5 mm, 45 cm, pistol grip
N3811630	SB-0535PC	SONICBEAT 5 mm, 35 cm, pistol grip
N3811730	SB-0545IC	SONICBEAT 5 mm, 45 cm, inline grip
N3811830	SB-0535IC	SONICBEAT 5 mm, 35 cm, inline grip
N3811930	SB-0520IC	SONICBEAT 5 mm, 20 cm, inline grip
N3812030	SB-0510IC	SONICBEAT 5 mm, 10 cm, inline grip



Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.