

ELECTROSURGICAL GENERATOR

ESG-400

Performance, Safety and Simplicity for All Your Surgical Disciplines



ESG-400 for Urological Applications


Intuitive User Interface

- Touch screen for easily finding and choosing your settings.
- Plug and play with instrument recognition via the UNIVERSAL socket.
- Automatic saline detection for ensuring you are using the right fluid.




Excellent and High Performance

- Tissue-adaptive power output for automatic application of optimal amount of energy.
- Advanced bipolar energy for continuous plasma activation with a 100% duty cycle during vaporization.
- High-powered cut support for reducing risk of undesired thermal spread.



Versatile Platform in the OR

- For open, lap and endoscopic surgery procedures.
- Fully equipped HF generator: for all surgical disciplines with monopolar, bipolar and advanced bipolar modes.



Specifications

Size and weight	370 × 465 × 156 mm, 12.5 kg
High frequency functions	Monopolar/bipolar
Max. high frequency power	320 W
Duty cycle	All modes: 25% RFCoag: 100%



Specifications, design, and accessories are subject to change without any notice or obligation on the part of the manufacturer.

OLYMPUS

OLYMPUS EUROPA SE & CO. KG
Postbox 10 49 08, 20034 Hamburg, Germany
Wendenstrasse 14-18, 20097 Hamburg, Germany
Phone: +49 40 23773-0, Fax: +49 40 233765
www.olympus-europa.com