


Clear View of the Surgical Site, Anytime

InstaClear: Lens Cleaner System



Clear View, Smooth Procedure

Safe and Precise Surgeries without Interruptions



Time Savings

- Combined suction and irrigation for better and faster lens cleaning.
- Clear view of the surgical site at all times.



Patient Safety

- Cleaning without scope removal.
- Design that reduces the duration of the procedure.



Flexibility

- Compatibility with all 4 mm HD and Ultra Olympus endoscopes as well as 4 mm Storz endoscopes.
- 0°, 30°, 45° and 70° working angles.
- Top or bottom light-cable positions.



Easy Handling


- Hands-free pedal activation.
- Adjustable flow rate for customized performance and surgical-site irrigation.

InstaClear: The Clear Choice

The InstaClear lens cleaner system cleans and clears the lens of a rigid endoscope almost instantly during surgery without having to remove it from the nasal cavity. The system includes a console that can be mounted on a standard IV pole for easy storage and transportation in the OR, a pedal that enables hands-free activation, a tube set that irrigates and cleans the lens, a suction function that clears the lens, and a surgical-field irrigation function. Additionally, the InstaClear comes with a series of easy-to-mount sheaths that accommodate 4 mm Olympus HD scopes, Olympus Ultra scopes and Storz 4 mm rigid endoscopes.


For more information, please visit:

➔ www.olympus.eu/instaclear



Also available for scopes with the light source on the top.

Accommodate All Viewing Angles and Light-Cable Positions



The InstaClear console can be used in combination with the DIEGO ELITE system or as a stand-alone lens-cleaning system.



Did You Know ...

Olympus Offers a Full See-and-Treat Portfolio for Office- and OR-Based Rhinology Procedures

Imaging

Treatment



E0492472EN · 11/19 · OEGK

Olympus reserves the right of errors, modification and changes of the service and/or product offerings.

OLYMPUS

OLYMPUS EUROPA SE & CO. KG
Postbox 10 49 08, 20034 Hamburg, Germany
Wendenstrasse 14-18, 20097 Hamburg, Germany
Phone: +49 40 23773-0, Fax: +49 40 233765
www.olympus-europa.com