

## moresolution Laparoscopic Power Morcellator

Tissue Extraction



**WARNING:** Uterine tissue may contain unsuspected cancer. The use of laparoscopic power morcellators during fibroid surgery may spread cancer, and decrease the long-term survival of patients. This information should be shared with patients when considering surgery with the use of these devices.

**WARNING:** Laparoscopic power morcellators should only be used with a technically compatible containment system.


Laparoscopic power morcellators are contraindicated in gynecologic surgery in which the tissue to be morcellated is known or suspected to contain malignancy. Laparoscopic power morcellators are contraindicated for removal of uterine tissue containing suspected fibroids in patients who are

- post-menopausal or over 50 years of age
- candidates for en bloc tissue removal through the vagina or via a mini-laparotomy

This content does not replace the instruction for use of the product.


# Tissue Extraction System

## Features and Benefits



### Single Use Cutting Modules and Valve Modules

Single use Cutting Modules for consistent cutting performance in every procedure



Single use Valve Module incl. Reducer Ø 10 mm to 5 mm for constant pressure even during instrument changes



### Dedicated Forceps

Available as Tenaculum and 2 x 3 Claw Forceps



330 mm or 420 mm working length



# moresolution Laparoscopic Power Morcellator

For more information, please visit

 **[www.olympus.eu/gynecology](http://www.olympus.eu/gynecology)**

moresolution® is a registered trademark of TROKAMED GmbH Germany.

The BLUE ENDO mark and logo are registered trademarks of blueendo GmbH international.

As medical knowledge is constantly growing, technical modifications or changes of the product design, product specifications and accessories may be required.

---

**OLYMPUS**

**OLYMPUS EUROPA SE & CO. KG**  
Postbox 10 49 08, 20034 Hamburg, Germany  
Wendenstrasse 20, 20097 Hamburg, Germany  
Phone: +49 40 23773-0, Fax: +49 40 233765  
[www.olympus-europa.com](http://www.olympus-europa.com)