

HANAROSTENT® Upper Esophageal Stent, Covered (CCC)							
	Model	Article Number	Stent (mm)		Total Length	Delivery Device	
			Working Length	Ø		Length (mm)	Ø (mm/Fr)
41188	ECP-18080-X070	E0422127	80	24-18-26	700	6/18	
	ECP-18080-X070	E0422128	52	80	24-18-26	700	6/18
	ECP-18100-X070	E0422129	72	100	24-18-26	700	6/18
	ECP-18120-X070	E0422130	92	120	24-18-26	700	6/18
	ECP-18150-X070	E0422131	122	150	24-18-26	700	6/18
	ECP-18170-X070	E0422132	142	170	24-18-26	700	6/18

HANAROSTENT® Upper Esophageal Stent, Covered (CCC)							
	Model	Article Number	Stent (mm)		Total Length	Delivery Device	
			Working Length	Ø		Length (mm)	Ø (mm/Fr)
41196	ECN-18-060-X070	E0422000	30	60	24-18-24	700	6/18
	ECN-18-090-X070	E0422001	60	90	24-18-24	700	6/18
	ECN-18-110-X070	E0422002	80	110	24-18-24	700	6/18

HANAROSTENT® Esophageal Stent, Covered (CCC)							
	Model	Article Number	Stent (mm)		Total Length	Delivery Device	
			Working Length	Ø		Length (mm)	Ø (mm/Fr)
41216	HES-18-060-070	E0422133	20	60	24-18-24	700	6/18
	HES-18-080-070	E0422134	40	80	24-18-24	700	6/18
	HES-18-110-070	E0422135	70	110	24-18-24	700	6/18
	HES-18-140-070	E0422136	100	140	24-18-24	700	6/18
	HES-18-170-070	E0422137	130	170	24-18-24	700	6/18
	HES-22-060-070	E0422138	20	60	28-22-28	700	8/24
	HES-22-080-070	E0422139	40	80	28-22-28	700	8/24
	HES-22-110-070	E0422140	70	110	28-22-28	700	8/24
	HES-22-140-070	E0422141	100	140	28-22-28	700	8/24
	HES-22-170-070	E0422142	130	170	28-22-28	700	8/24
	HES-22-180-070	E0422143	140	180	28-22-28	700	8/24

HANAROSTENT® Esophageal Stent, Covered (CCC)							
	Model	Article Number	Stent (mm)		Total Length	Delivery Device	
			Working Length	Ø		Length (mm)	Ø (mm/Fr)
41200	ECV-20-080-X070	E0422003	40	80	26-20-26	700	6/18
	ECV-20-110-X070	E0422004	70	110	26-20-26	700	6/18
	ECV-20-140-X070	E0422005	100	140	26-20-26	700	6/18
	ECV-20-170-X070	E0422006	130	170	26-20-26	700	6/18

HANAROSTENT® Esophageal Stent, Covered (CCC)							
	Model	Article Number	Stent (mm)		Total Length	Delivery Device	
			Working Length	Ø		Length (mm)	Ø (mm/Fr)
41232	NES-28-060-070	E0422144	20	60	34-28-34	700	8/24
	NES-28-080-070	E0422145	40	80	34-28-34	700	8/24
	NES-28-100-070	E0422146	60	100	34-28-34	700	8/24
	NES-28-110-070	E0422007	70	110	34-28-34	700	8/24
	NES-28-120-070	E0422147	80	120	34-28-34	700	8/24
	NES-28-140-070	E0422008	100	140	34-28-34	700	8/24
	NES-28-170-070	E0422009	130	170	34-28-34	700	8/24

HANAROSTENT® Esophageal Stent, Partially Covered (NCN)							
	Model	Article Number	Stent (mm)		Total Length	Delivery Device	
			Working Length	Ø		Length (mm)	Ø (mm/Fr)
41204	EPB-18060-X070	E0422148	20	60	24-18-24	700	6/18
	EPB-18080-X070	E0422149	40	80	24-18-24	700	6/18
	EPB-18110-X070	E0422150	70	110	24-18-24	700	6/18
	EPB-18140-X070	E0422151	100	140	24-18-24	700	6/18
	EPB-18170-X070	E0422152	130	170	24-18-24	700	6/18
	EPB-22060-Z070	E0422155	20	60	28-22-28	700	8/24
	EPB-22080-Z070	E0422156	40	80	28-22-28	700	8/24
	EPB-22090-Z070	E0422157	50	90	28-22-28	700	8/24
	EPB-22110-Z070	E0422158	70	110	28-22-28	700	8/24
	EPB-22120-Z070	E0422159	80	120	28-22-28	700	8/24
	EPB-22140-Z070	E0422160	100	140	28-22-28	700	8/24
	EPB-22160-Z070	E0422161	120	160	28-22-28	700	8/24
	EPB-22170-Z070	E0422162	130	170	28-22-28	700	8/24

HANAROSTENT® Esophageal Stent, Partially Covered (NCN)							
	Model	Article Number	Stent (mm)		Total Length	Delivery Device	
			Working Length	Ø		Length (mm)	Ø (mm/Fr)
41212	EPW-20-080-070	E0422010	40	80	26-20-26	700	6/18
	EPW-20-100-070	E0421072	60	100	26-20-26	700	6/18
	EPW-20-110-070	E0422011	70	110	26-20-26	700	6/18
	EPW-20-140-070	E0422012	100	140	26-20-26	700	6/18
	EPW-20-150-070	E0421073	110	150	26-20-26	700	6/18

HANAROSTENT® Esophageal Stent, Skidproof, Covered (CCC)							
	Model	Article Number	Stent (mm)		Total Length	Delivery Device	
			Working Length	Ø		Length (mm)	Ø (mm/Fr)
41194	ECD-20-080-070	E0422169	40	80	26-20-26	700	6/18
	ECD-20-120-070	E0422170	80	120	26-20-26	700	6/18
	ECD-20-140-070	E0422171	100	140	26-20-26	700	6/18
	ECD-20-180-070	E0422172	140	180	26-20-26	700	8/24


HANAROSTENT® Esophageal Stent, Uncovered (NNN)							
	Model	Article Number	Stent (mm)		Total Length	Delivery Device	
			Working Length	Ø		Length (mm)	Ø (mm/Fr)
41220	HNS-22-090-070	E0422173	54	90	28-22-28	700	4/12
	HNS-22-120-070	E0422174	84	120	28-22-28	700	4/12
	HNS-22-160-070	E0422175	124	160	28-22-28	700	4/12

HANAROSTENT® Esophageal Stent, Covered, for Benign Strictures (CCC)							
	Model	Article Number	Stent (mm)		Total Length	Delivery Device	
			Working Length	Ø		Length (mm)	Ø (mm/Fr)
41186	EBN22060-Z070	E0422191	20	60	28-22-28	700	8/24
	EBN22080-Z070	E0422015	40	80	28-22-28	700	8/24
	EBN22100-Z070	E0422192	60	100	28-22-28	700	8/24
	EBN22120-Z070	E0422016	80	120	28-22-28	700	8/24

HANAROSTENT® Esophageal Stent, Covered (CCC), TTS							
	Model	Article Number	Stent (mm)		Total Length	Delivery Device	
			Working Length	Ø		Length (mm)	Ø (mm/Fr)
86560	ECBA-18-060-230	E0420533	30	60	26-18-26	2300	3.5/10.5
	ECBA-18-090-230	E0420534	60	90	26-18-26	2300	3.5/10.5
	ECBA-18-100-230	E0420535	70	100	26-18-26	2300	3.5/10.5
	ECBA-18-120-230	E0420536	90	120	26-18-26	2300	3.5/10.5
	ECBA-18-150-230	E0420537	120	150	26-18-26	2300	3.5/10.5

HANAROSTENT® Esophageal Stent, Partially Covered (NCN), TTS							
	Model	Article Number	Stent (mm)		Total Length	Delivery Device	
			Working Length	Ø		Length (mm)	Ø (mm/Fr)
41206	EPBA-20-060-230	E0422198	30	60	26-20-26	2300	3.5/10.5
	EPBA-20-080-230	E0422199	50	80	26-20-26	2300	3.5/10.5
	EPBA-20-090-230	E0422200	60	90	26-20-26	2300	3.5/10.5
	EPBA-20-100-230	E0422200	70	100	26-20-26	2300	3.5/10.5
	EPBA-20-120-230	E0422201	90	120	26-20-26	2300	3.5/10.5
	EPBA-20-150-230	E0422202	120	150	26-20-26	2300	3.5/10.5

HANAROSTENT® Esophageal Stent with Anti-Reflex Function, Covered (CCC)							
	Model	Article Number	Stent (mm)		Total Length	Delivery Device	
			Working Length	Ø		Length (mm)	Ø (mm/Fr)
41218	HESV-18-090-070	E0422177	90	120	24-18-24	700	8/24
	HESV-18-120-070	E0422177	80	120	24-18-24	700	8/24
	HESV-18-160-070	E0422178	120	160	24-18-24	700	8/24
	HESV-22-060-070	E0422179	20	60	28-22-28	700	8/24
	HESV-22-090-070	E0422180	50	90	28-22-28	700	8/24
	HESV-22-120-070	E0422181	80	120	28-22-28	700	8/24
	HESV-22-160-070	E0422182	120	160	28-22-28	700	8/24
	HESV-22-170-070	E0422183	130	170	28-22-28	700	8/24

HANAROSTENT® Esophageal Stent with Anti-Reflex Function, Covered (CCC)							
	Model	Article Number	Stent (mm)		Total Length	Delivery Device	
			Working Length	Ø		Length (mm)	Ø (mm/Fr)
41202	ECW-22-090-Z070	E0422013	50	90	28-22-28	700	8/24
	ECW-22-120-Z070	E0422014	80	120	28-22-28	700	8/24
	ECW-22-160-Z070	E0420582	120	160	28-22-28	700	8/24