


Available in two sizes to fit all types of endoscope.

Protects the endoscope during transportation and storage to reduce the risk of damage and increase equipment availability.

Minimal contact with the endoscope, allowing aeration and facilitating easy removal.

Unique clean/dirty tab helps to prevent reuse with subsequent endoscopes.

With sterile packaging and only suitable for single use to reduce cross contamination.



PROTECH

Flexible Endoscope Tip Protector

A simple, safe and highly effective method of protecting the delicate optics in the distal end of the endoscope.

Ordering Information

Article Number	Product Description	Unit
E0425694	PROTECH S	100 pieces/box
E0425695	PROTECH L	100 pieces/box
E0420577	PROTECH XL	50 pieces/box



DISTRIBUTOR

OLYMPUS EUROPA SE & CO. KG
 Postbox 10 49 08, 20034 Hamburg, Germany
 Wendenstrasse 20, 20097 Hamburg, Germany
 Phone: +49 40 23773-0, Fax: +49 40 233765
 www.olympus-europa.com

MANUFACTURER

Meditech Endoscopy Limited, Bowden House,
 Luckyn Lane, Basildon, Essex, SS14 3AX, UK
 Phone: +44 1268 534035
 Email: admin@i-medical.net
 www.meditechendoscopy.com