

OLYMPUS

Your Vision, Our Future

LIGHT GUIDE CABLES

Eco-friendly and durable



OPTIMAL LIGHT TRANSMISSION

Ensuring you have a good view can be the key to success – and this particularly applies to endoscopic procedures. Olympus light guide cables mean you can focus on your day-to-day work, and not have to worry about light being correctly transmitted.

When combined with an Olympus light source and the required scope, you'll benefit from optimal illumination to suit your needs.


- Eco-friendly in terms of RoHS compliance
- CF applied part classification for use in cardiac surgery
- Optimal light transmission thanks to fused optical fibers
- Durability: optimized tube and fused fiber bundles for extremely high durability against mechanical stress and autoclaving
- High temperature resistance
- Compatible with a wide variety of endoscopes if used with suitable adapters
- Autoclavable



Specifications


Delivery of all light-guide cables includes light-guide adapter WA03320A.

WA03300A Light-guide cable, 2.8 mm, 3 m, CF type



maximum insertion portion width of the endoscope ≤ 4.1 mm

WA03310A Light-guide cable, 4.25 mm, 3 m, CF type



maximum insertion portion width of the endoscope > 4.1 mm

Spare Parts

WA03320A Spare light-guide adapter, Olympus light-guide cable to Olympus light source

