

OLYMPUS

Your Vision, Our Future

COMPACT TROLLEY TC-C4



TC-C4 COMPACT TROLLEY

The TC-C4 Compact Trolley is designed as a compact trolley but which is large enough to be used in conjunction with EVIS Exera III imaging system.

Compact

The compact design means the TC-C4 is ideal for use in the outpatient or doctor's office where space is limited.

Easy to Move

The compact trolley is easy to move due to its light weight and integrated handles at the front and rear of the top tray.

Monitor Support

The TC-C4 is offered with a choice of monitor arm or monitor mount allowing optimised procedure image viewing position.

Large Top Tray

The top tray provides a large space for working or can be used to place equipment to increase storage capacity.

Cable management

Optimised cable management is achieved through channels in the columns and cable ties on rear of shelf.

Integrated Separation Transformer

Optional transformer placed under the base meaning increased shelf space for equipment loading.

Keyboard tray

A sliding keyboard tray is located under the top shelf.



ACCESSORIES

MAJ-2146

Sliding Keyboard Tray



MAJ-1684

Scope Pole Kit



MAJ-2148

Rear Panel



MAJ-1686

Drawer Unit



MAJ-2147

Compact LCD Monitor Arm



MAJ-1683

LCD Monitor Mount



MAJ-1656

Footswitch Holder



MAJ-1639

CO₂ Cylinder Holder (140 mm)



MAJ-1665

Camera Head Holder



MAJ-1663

EUS Arm Mount Kit



MAJ-1658

Universal Storage Container



MAJ-110

Jar Mounting Kit TC



TC-C4 STANDARD SETS

Standard Sets include:

TC-C4 Compact Trolley plus items marked (●)

Accessories available as an optional extra (o)

MAJ Description	Part Description	Part Number	Region	Set 1	Set 2	Set 3	Set 4
			EUROPE	K10029882	K10029887	K10029892	K10029897
			JAPAN	K10029885	K10029890	K10029895	K10029900
			UK	K10029881	K10029886	K10029891	K10029896
			US	K10029884	K10029889	K10029894	K10029899
			ROW	K10029883	K10029888	K10029893	K10029898
MAJ-2146	Sliding Keyboard Tray	K10027573	O	●	●	●	●
MAJ-1684	Scope Pole Kit	K10027281	O	●	●		
MAJ-2148	Rear Panel	K10027574	O	●	●	●	
	Transformer			●		●	
MAJ-1686	Drawer Unit	K10027569	O				
MAJ-2147	Compact LCD Monitor Arm	K10026879	O				
MAJ-1683	LCD Monitor Mount	K10026888	O				
MAJ-1656	Footswitch Holder	K10021790	O				
MAJ-1639	CO ₂ Cylinder Holder (140 mm)	K10021041	O				
MAJ-1665	Camera Head Holder	K10022056	O				
MAJ-1663	EUS Arm Mount Kit	K10021797	O				
MAJ-1658	Universal Storage Container	K10021792	O				
MAJ-110	Jar Mounting Kit TC	K7503556	O				
MAJ-1654	Equipotential Terminal Strip	K10021352	O				
	Blister Pack Toolkit		O	●	●	●	●
	IEC Lead Set		O	●		●	

Specifications

TC-C4

Dimensions	1110 mm (H) x 675 mm (D) x 534 mm (W)
	Height to top tray work surface - 1100 mm
	Height with 26" LCD monitor in highest position - 1945 mm
Weight	Maximum 87 kg (unladen) including fitted transformer

Voltage

The TC-C4 can be supplied with a separating transformer based on one of the following voltage options:

MAJ-1646	WM-T2 Transformer 100V
MAJ-1647	WM-T2 Transformer 110-120V
MAJ-1648	WM-T2 Transformer 220-240V

Construction

Trolley

Columns	Extruded profiled aluminium
Shelves	Mild steel
Top Tray	Injection moulded ABS

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

OLYMPUS

OLYMPUS KEYMED
KeyMed House, Stock Road, Southend-on-Sea,
Essex SS2 5QH, United Kingdom
Telephone: +44 (0)1702 616333 Fax: +44 (0)1702 465677
e-mail: keymed@keymed.co.uk www.keymed.co.uk

For a complete listing of
sales and distribution locations visit:
www.olympus.com