


Irrigation and Insufflation Tubing and Adaptors

Selection Guide



Irrigation and Insufflation Tubing and Adaptors

Selection Guide



* Available from other vendors

Irrigation and Insufflation Tubing and Adaptors

Selection Guide

Combined
Irrigation
&
Insufflation
Tube set
(Hybrid)

Irrigation Tubing with Air (includes lens cleaning)		
Article Number	Article Name	Compatible Water Bottles (1000 ml)
K10035004		Baxter® UKF7114 / AKF7114
		Baxter® 2F7114
		Baxter® AHF7114
		B. Braun® R5000-01
		Hospira® 07139-09
K10035005		Baxter® UKF7114 / AKF7114
		B. Braun Ecotainer® 0082479E
		Fresenius Kabi® 3158589 / B230532
		Fresenius Kabi Ampuwa® 1088813
		Laboratoire Aguettant Versol® 600500


Irrigation Tubing with CO ₂ (includes lens cleaning)		
Article Number	Article Name	Compatible Water Bottles (1000 ml)
K10035002		Baxter® UKF7114 / AKF7114
		Baxter® 2F7114
		Baxter® AHF7114
		B. Braun® R5000-01
		Hospira® 07139-09
K10035003		Baxter® UKF7114 / AKF7114
		B. Braun Ecotainer® 0082479E
		Fresenius Kabi® 3158589 / B230532
		Fresenius Kabi Ampuwa® 1088813
		Laboratoire Aguettant Versol® 600500



+



See System 7
(page 9)



+



See System 7
(page 9)



+




See System 8
(page 9)




+



See System 8
(page 9)




	CO ₂ Source Tube (Single day use)
	MAJ-1741 High Flow Gas Tube
	MAJ-1742 Low Flow Gas Tube
	MAJ-1816 Extra Low Flow Gas Tube
	CO ₂ Source Tubing Adaptor (Single day use)



Tube Holders
MAJ-2160
(WM-NP2, WM-DP2, WM-NP3, WM-DP3, TC-C3 & TC-C4 compatible)

Irrigation and Insufflation Tubing and Adaptors

Selection Guide




Sterile Water Holder
MAJ-2166


Irrigation and Insufflation Tubing and Adaptors

System Schematic Diagrams


System 1
Fully Disposable Irrigation Only (Auxiliary channel)




System 2
Semi-disposable Irrigation Only (Auxiliary channel)




System 5
Air Insufflation Only (Reusable)




System 6
CO₂ Insufflation Only (Reusable)




System 3
ESD Injection / Irrigation




System 4
Fully Disposable Irrigation Only (Instrument channel)



System 7
Combined Irrigation and Insufflation Tubeset with Air



System 8
Combined Irrigation and Insufflation Tubeset with CO₂




* CO₂ source can be positioned on shelf above Flushing Pump without the need for CO₂ Source Tube


Irrigation and Insufflation Tubing and Adaptors

System Schematic Diagrams

System 9
Air Insufflation Only



System 10
CO₂ Insufflation Only



Irrigation and Insufflation Tubing and Adaptors

Ordering Information				
Article Number	Article Name	Description	Note	Quantity Unit/Box
K10035002	MAJ-2207	Irrigation Tubing with CO ₂ Cap A	Single day	10
K10035003	MAJ-2208	Irrigation Tubing with CO ₂ Cap B	Single day	10
K10035004	MAJ-2209	Irrigation Tubing with Air Cap A	Single day	10
K10035005	MAJ-2210	Irrigation Tubing with Air Cap B	Single day	10
K10023086	MAJ-1651	Auxiliary Channel Water Tube Set	Single day	10
K10020736	MAJ-1652	Auxiliary Channel Adaptor	Single patient	100
K10001147	MAJ-1608	Auxiliary Channel Water Tube	Single day	50
K10016135	MAJ-1608	Auxiliary Channel Water Tube	Single day	10
N3499350	MAJ-855	Auxiliary Water Tube	Single patient	1
K10007072	MAJ-1606	Instrument Channel Adaptor	Single patient	100
K10016091	MAJ-1606	Instrument Channel Adaptor	Single patient	10
K10001146	MAJ-1607	Instrument Channel Water Tube	Single patient	50
K10016136	MAJ-1607	Instrument Channel Water Tube	Single patient	10
K10026641	MAJ-1681	Accessory Port Tube with Bottle Cap (Use with ESD Knives)	Single patient	10
K10026686	MAJ-1682	Accessory Port Tube with Saline Spike (Use with ESD Knives)	Single patient	10
N2654500	MAJ-1741	UCR High Flow Gas Tube	Reusable	1
N2654600	MAJ-1742	UCR Low Flow Gas Tube	Reusable	1
N3493330	MAJ-1816	UCR Extra Low Flow Gas Tube	Reusable	1
028595	MAJ-901	Insufflation Water Container (Air)	Reusable	1
028596	MAJ-902	Insufflation Water Container (CO ₂)	Reusable	1
K10028141	MAJ-2160	Irrigation Tube Holder (WM-P2, WM-DP2, TC-C3 & TC-C4 compatible)	Reusable (low level disinfection)	1
E0420575	Contu CO ₂	CO ₂ Connection Tube including removable UCR adapter	Single day	10

Irrigation and Insufflation Tubing and Adaptors

Ordering Information (continued)

Article Number	Article Name	Description	Note	Quantity Unit/Box
K10037161	MAJ-2220	Insufflation Tubing CO ²	Single day	10
K10037162	MAJ-2221	Insufflation Tubing CO ²	Single day	10
K10037163	MAJ-2222	Insufflation Tubing Air	Single day	10
K10037164	MAJ-2223	Insufflation Tubing Air	Single day	10
K10037404	–	ADA-CO ² -OUS - CO ² Source Tubing Adaptor	Single day	100
K10030036	MAJ-2166	Sterile Water Holder	Reusable (low level disinfection)	1

As medical knowledge is constantly growing, technical modifications or changes of the product design, product specifications, accessories and service offerings may be required.



OLYMPUS EUROPA SE & CO. KG
Postbox 10 49 08, 20034 Hamburg, Germany
Wendenstrasse 20, 20097 Hamburg, Germany
Phone: +49 40 23773-0, Fax: +49 40 233765
www.olympus-europa.com

For a complete listing of
sales and distribution locations visit:
www.olympus.com