

Esophagus

HANAROSTENT® ECN – Single-Use Upper Esophageal Covered Stents (CCC)						
	Model	Stent (mm)			Delivery Device	
		Working Length	Total Length	Ø	Length (mm)	Ø (mm/Fr)
41196	ECN-18-060-X070	30	60	24-18-24	700	6/18
	ECN-18-090-X070	60	90	24-18-24	700	6/18
	ECN-18-110-X070	80	110	24-18-24	700	6/18

HANAROSTENT® ECV – Single-Use Esophageal Covered Stents (CCC)						
	Model	Stent (mm)			Delivery Device	
		Working Length	Total Length	Ø	Length (mm)	Ø (mm/Fr)
41200	ECV-20-080-X070	40	80	26-20-26	700	6/18
	ECV-20-110-X070	70	110	26-20-26	700	6/18
	ECV-20-140-X070	100	140	26-20-26	700	6/18
	ECV-20-170-X070	130	170	26-20-26	700	6/18

HANAROSTENT® NES26 – Single-Use Esophageal Covered Stents (CCC)						
	Model	Stent (mm)			Delivery Device	
		Working Length	Total Length	Ø	Length (mm)	Ø (mm/Fr)
41232	NES-28-110-070	70	110	34-28-34	700	8/24
	NES-28-140-070	100	140	34-28-34	700	8/24
	NES-28-170-070	130	170	34-28-34	700	8/24

HANAROSTENT® EPW – Single-Use Esophageal Partially Covered Stents (NCN)						
	Model	Stent (mm)			Delivery Device	
		Working Length	Total Length	Ø	Length (mm)	Ø (mm/Fr)
41212	EPW-20-080-X070	40	80	26-20-26	700	6/18
	EPW-20-110-X070	70	110	26-20-26	700	6/18
	EPW-20-140-X070	100	140	26-20-26	700	6/18

HANAROSTENT® EBN – Single-Use Esophageal Covered Stents for Benign Strictures (CCC)						
	Model	Stent (mm)			Delivery Device	
		Working Length	Total Length	Ø	Length (mm)	Ø (mm/Fr)
41186	EBN-22-060-Z070	20	60	28-22-28	700	8/24
	EBN-22-080-Z070	40	80	28-22-28	700	8/24
	EBN-22-100-Z070	60	100	28-22-28	700	8/24
	EBN-22-120-Z070	80	120	28-22-28	700	8/24

HANAROSTENT® ECD – Single-Use Esophageal Covered, Skidproof Stents (CCC)						
	Model	Stent (mm)			Delivery Device	
		Working Length	Total Length	Ø	Length (mm)	Ø (mm/Fr)
41194	ECD-20-080-070	40	80	26-20-26	700	6/18
	ECD-20-120-070	80	120	26-20-26	700	6/18
	ECD-20-140-070	100	140	26-20-26	700	6/18
	ECD-20-180-070	140	180	26-20-26	700	6/18

HANAROSTENT® HNS – Single-Use Esophageal Uncovered Stents (NNN)						
	Model	Stent (mm)			Delivery Device	
		Working Length	Total Length	Ø	Length (mm)	Ø (mm/Fr)
41220	HNS22-090-070	50	90	28-22-28	700	4/12
	HNS22-120-070	80	120	28-22-28	700	4/12
	HNS22-160-070	120	160	28-22-28	700	4/12

HANAROSTENT® EPBA – Single-Use Esophageal Partially Covered Stents (NCN), TTS						
	Model	Stent (mm)			Delivery Device	
		Working Length	Total Length	Ø	Length (mm)	Ø (mm/Fr)
41206	EPBA-18-060-230	30	60	26-18-26	2300	3.4/10.5
	EPBA-18-080-230	50	80	26-18-26	2300	3.4/10.5
	EPBA-18-100-230	70	100	26-18-26	2300	3.4/10.5
	EPBA-18-120-230	90	120	26-18-26	2300	3.4/10.5
	EPBA-18-150-230	120	150	26-18-26	2300	3.4/10.5
	EPBA-20-060-230	30	60	26-20-26	2300	3.4/10.5
	EPBA-20-080-230	60	80	26-20-26	2300	3.4/10.5
	EPBA-20-090-230	60	90	26-20-26	2300	3.4/10.5
	EPBA-20-120-230	90	120	26-20-26	2300	3.4/10.5
	EPBA-20-150-230	120	150	26-20-26	2300	3.4/10.5

HANAROSTENT® ECW – Single-Use Esophageal Covered Stents with Anti-Reflux Function (CCC)						
	Model	Stent (mm)			Delivery Device	
		Working Length	Total Length	Ø	Length (mm)	Ø (mm/Fr)
41202	ECW-22-090-Z070	50	90	28-22-28	700	8/24
	ECW-22-120-Z070	80	120	28-22-28	700	8/24

HANAROSTENT® ECBS – Single-Use Esophageal Covered Stents (CCC)						
	Model	Stent (mm)			Delivery Device	
		Working Length	Total Length	Ø	Length (mm)	Ø (mm/Fr)
41192	ECBS-26-180-090	130	180	34-26-36	900	8/24
	ECBS-26-210-090	160	210	34-26-36	900	8/24
	ECBS-26-240-090	190	240	34-26-36	900	8/24

HANAROSTENT® ECBB – Single-Use Esophageal Covered Stents (CCC)						
	Model	Stent (mm)			Delivery Device	
		Working Length	Total Length	Ø	Length (mm)	Ø (mm/Fr)
41188	ECBB-28-180-090	145	180	32-28-28	900	8/24
	ECBB-28-210-090	175	210	32-28-28	900	8/24
	ECBB-28-240-090	205	240	32-28-28	900	8/24

Biliary System

HANAROSTENT® BNHS – Single-Use Hilar Biliary Uncovered Stents (NNN), TTS						
	Model	Stent (mm)			Delivery Device	
		Working Length	Total Length	Ø	Length (mm)	Ø (mm/Fr)
41158	BNHS-06-040-180	40	40	6	1800	1.9/7.5
	BNHS-06-060-180	60	60	6	1800	1.9/7.5
	BNHS-06-080-180	80	80	6	1800	1.9/7.5
	BNHS-06-100-180	100	100	6	1800	1.9/7.5
	BNHS-06-120-180	120	120	6	1800	1.9/7.5
	BNHS-08-040-180	40	40	8	1800	1.9/7.5
	BNHS-08-060-180	60	60	8	1800	1.9/7.5
	BNHS-08-080-180	80	80	8	1800	1.9/7.5
	BNHS-08-100-180	100	100	8	1800	1.9/7.5
	BNHS-08-120-180	120	120	8	1800	1.9/7.5


HANAROSTENT® BNL – Single-Use Hilar Biliary Uncovered Stents (NNN), TTS						
	Model	Stent (mm)			Delivery Device	
		Working Length	Total Length	Ø	Length (mm)	Ø (mm/Fr)
41160	BNL-10-060-180	46	60	13.5-10-13.5	1800	2.33/7
	BNL-10-080-180	66	80	13.5-10-13.5	1800	2.33/7

HANAROSTENT® BNY – Single-Use Hilar Biliary Uncovered Stents (NNN), TTS						
	Model	Stent (mm)			Delivery Device	
		Working Length	Total Length	Ø	Length (mm)	Ø (mm/Fr)
41162	BNY-10-060-180	46	60	13.5-10-13.5	1800	2.33/7
	BNY-10-080-180	66	80	13.5-10-13.5	1800	2.33/7
	BNY-10-100-180	86	100	13.5-10-13.5	1800	2.33/7
	BNY-10-120-180	106	120	13.5-10-13.5	1800	2.33/7

HANAROSTENT® BNY – Single-Use Hilar Biliary Uncovered Stents for Percutaneous Application (NNN)						
	Model	Stent (mm)			Delivery Device	
		Working Length	Total Length	Ø	Length (mm)	Ø (mm/Fr)
41162	BNY-10-060-060	46	60	13.5-10-13.5	600	2.33/7
	BNY-10-080-060	66	80	13.5-10-13.5	600	2.33/7
	BNY-10-100-060	86	100	13.5-10-13.5	600	2.33/7
	BNY-10-120-060	106	120	13.5-10-13.5	600	2.33/7

HANAROSTENT® SHS – Single-Use Biliary Uncovered Stents (NNN), TTS						
	Model	Stent (mm)			Delivery Device	
		Working Length	Total Length	Ø	Length (mm)	Ø (mm/Fr)
41238	SHS-10-040-180	26	40	13.5-10-13.5	1800	2.33/7
	SHS-10-060-180	46	60	13.5-10-13.5	1800	2.33/7
	SHS-10-080-180	66	80	13.5-10-13.5	1800	2.33/7
	SHS-10-100-180	86	100	13.5-10-13.5	1800	2.33/7

HANAROSTENT® SHS – Single-Use Biliary Uncovered Stents (NNN), TTS – Distal Wire-Guided, FAST™ Delivery System						
	Model	Stent (mm)			Delivery Device	
		Working Length	Total Length	Ø	Length (mm)	Ø (mm/Fr)
41238	SHS-10-040-S180	26	40	13.5-10-13.5	1800	2.83/8.5
	SHS-10-060-S180	46	60	13.5-10-13.5	1800	2.83/8.5
	SHS-10-080-S180	66	80	13.5-10-13.5	1800	2.83/8.5
	SHS-10-100-S180	86	100	13.5-10-13.5	1800	2.83/8.5



HANAROSTENT® SHS – Single-Use Hilar Biliary Uncovered Stents for Percutaneous Application (NNN)						
	Model	Stent (mm)			Delivery Device	
		Working Length	Total Length	Ø	Length (mm)	Ø (mm/Fr)
41238	SHS-10-040-060	26	40	13.5-10-13.5	600	2.33/7
	SHS-10-060-060	46	60	13.5-10-13.5	600	2.33/7
	SHS-10-080-060	66	80	13.5-10-13.5	600	2.33/7
	SHS-10-100-060	86	100	13.5-10-13.5	600	2.33/7

HANAROSTENT® BNA – Single-Use Biliary Uncovered Stents (NNN), TTS						
	Model	Stent (mm)			Delivery Device	
		Working Length	Total Length	Ø	Length (mm)	Ø (mm/Fr)
 41160	BNA-08-040-G180	26	40	11.5-8-11.5	1800	2.66/8
	BNA-08-050-G180	36	50	11.5-8-11.5	1800	2.66/8
	BNA-08-060-G180	46	60	11.5-8-11.5	1800	2.66/8
	BNA-08-080-G180	66	80	11.5-8-11.5	1800	2.66/8
	BNA-08-100-G180	86	100	11.5-8-11.5	1800	2.66/8
	BNA-10-040-G180	26	40	13.5-10-13.5	1800	2.66/8
	BNA-10-050-G180	36	50	13.5-10-13.5	1800	2.66/8
	BNA-10-060-G180	46	60	13.5-10-13.5	1800	2.66/8
	BNA-10-080-G180	66	80	13.5-10-13.5	1800	2.66/8
	BNA-10-100-G180	86	100	13.5-10-13.5	1800	2.66/8