

# Hytrack 2


Managing Reprocessing from a New Perspective



# Hytrack

## Future Proof Software Solution for an Efficient Reprocessing Workflow

Hytrack manages your reprocessing from a new perspective by enabling you to track the reprocessing status of your endoscopes at any time and perform and trace the correct reprocessing to meet regulatory or hospital specific requirements.



## Future Proof Software Solution for an Efficient Reprocessing Workflow



### Full Traceability and Transparency

Provides a comprehensive overview of the current status of your machines and endoscopes, so you can see your endoscope's reprocessing stage and availability for the next procedure at any time.



### Protecting User Privacy

Places a high priority on data security and meeting GDPR requirements.



### Ease of Use

With an intuitive touch enabled and app based user interface, Hytrack supports you in documenting your reprocessing process step-by-step. With fewer interactions thanks to the RFID system, it offers peace of mind for your staff.



### Efficient Management

Minimizes errors and accelerates administration of hospital and loaner endoscopes. Save administration time with the seamless integration of Olympus devices.



### Compliance with Regulations

Manifold configuration options allow you to meet tailored regulatory or hospital specific hygiene documentation requirements to help you to adhere to endoscope reprocessing regulations.



### Maximized Uptime

Specific and customizable reports and an integrated task scheduling function enable you to optimize the performance of the reprocessing department.

# Hytrack

## Versatile Modularity Meeting Your Needs

### License Options for Hytrack Modules

Hytrack supports the reprocessing process by keeping track of all reprocessing steps.

#### Basic License

##### Foundation of Your Reprocessing Documentation

Documents and keeps track of all automated reprocessing steps by linking your endoscope pool with all connected devices and users. A system interface to Endobase or other non-Olympus documentation systems allows seamless integration into your hospital IT environment.



## Versatile Modularity Meeting Your Needs

The modular structure through licenses building on the Basic License allows you to select the options that best meet the needs of your department. Become even more flexible by subscribing the **additional licenses**:

### Documentation License

Prevent faulty reprocessing by creating your own or regulatory required procedural instructions for reprocessing documentation, document reprocessing with just a few clicks and having the documentation at your fingertips.

### Endoscope Dashboard License

Stop spending time on searching for endoscopes, you can see which reprocessing stage your endoscope is in and how long it will take until the documentation will be completed.

### Quick Scan License

Less stress at reprocessing stations, including the manual ones, thanks to faster and more convenient documentation enabled by RFID technology – all without touching a screen, mouse or keyboard.

### Machine Dashboard License

Optimize machine utilization by always having an overview of the current status and load of your reprocessing devices.

### Reports License

Bring the performance of your reprocessing department to the next level with specific and customizable reports.

### Procedure Interface License

Stay on top of your upcoming procedure by reconciling whether your endoscopes are available or not – thanks to the handy interface to patient scheduling systems.

# Hytrack

## Meets the IT Standards of a Modern Hospital

- ✓ Hytrack features low system requirements and supports current standards for client and server hardware and software\*.
- ✓ Hytrack supports the obligation of the operator to prove hospital IT cybersecurity, e.g., through background logging of system, machine, and user data stored and protected on the hospital network.
- ✓ Existing user directories of the hospital can be used consistently thanks to the Hytrack LDAP interface. Hytrack is designed with an API interface to be integrated into other non-Olympus software.
- ✓ Hytrack can be administered centrally to reduce the workload for the IT department.

\*Please contact your Olympus representative for further information.



**Complementary Service Option for Hytrack and your reprocessing equipment.**

Clinicians need to be able to rely on their equipment to provide high-quality patient care every day. By choosing Olympus services, customers can be assured that they will receive the most value over the contract term and beyond product guarantees. Olympus offers needs-based modular service contracts to suit the hospital's exact requirements.



**i** Your Olympus partner will be happy to assist you in finding the right service solution for your environment.



**Assure Optimal Equipment Availability and Budget Security**

**Equipment Availability** 

On site or remote interventions enable reliable procedure scheduling and revenue forecasting.

**Budget Security** 

Full cost transparency and reduced administration effort thanks to fixed annual service budgets.

**Patient and User Safety** 

Latest service information, specialized tooling, and testing methods to ensure equipment remains safe for you and your patients.

As medical knowledge is constantly growing, technical modifications or changes of the product design, product specifications and accessories may be required.