


STONE MANAGEMENT

# No-Tip Nitinol Baskets

Flex-Catch, Ultra-Catch and X-Catch for Reliable Stone Removal



# No-Tip Nitinol Baskets

## Versatile Portfolio

### Flex-Catch 1.9 Fr

The Flex-Catch NT stone retrieval device is flexible enough to conform to the anatomy while maintaining its shape to facilitate capture and release of stones. The no-tip design with four round wires minimizes tissue trauma during manipulation. This new basket is designed such that it gradually opens to capture stones of varying sizes.

- **Small 1.9 Fr. sheath:** Improved irrigation and visibility within small working channels.
- **Flexible composite sheath:** Maximum scope deflection and smooth scope passage.
- **Gradual basket opening:** Easy control over basket size to capture even small fragments.


### Ultra-Catch 4-Wire No-Tip Nitinol Basket

- Torque-responsive shaft and ergonomic handle provide the control and security to ensure safe and effective stone removal.
- Rotation wheel of the handle allows for single-handed operation.
- Highly durable twisted-wire cage design maintains its shape and form even under extreme duress.
- Nitinol construction makes the shaft kink-resistant, thereby delivering safety, functionality and durability.
- Flexible distal end minimizes deflection loss of the scope.


### X-Catch 6-Wire No-Tip Nitinol Basket

- Paired wire basket for higher radial dilation force for stone manipulation in restricted anatomy.
- Paired construction enables capture of several small stones in one pass to save time.
- No-tip design for atraumatic stone capture in the kidney.

## Flex-Catch Features



## Ultra-Catch and X-Catch Features



# No-Tip Nitinol Baskets

 [www.olympus.eu/medical](http://www.olympus.eu/medical)



E0492260 · 500 · 09/18 · OELG · ED

Specifications, design, and accessories are subject to change without any notice or obligation on the part of the manufacturer.

**OLYMPUS**

**OLYMPUS EUROPA SE & CO. KG**

Postbox 10 49 08, 20034 Hamburg, Germany  
Wendenstrasse 14-18, 20097 Hamburg, Germany  
Phone: +49 40 23773-0, Fax: +49 40 233765  
[www.olympus-europa.com](http://www.olympus-europa.com)